

↓ 400-901-9800⋈ sales@bioss.com.cn⋈ techsupport@bioss.com.cn

## Mouse Anti-OCT2 (ready to ues)

Cat. Number: BH0136

Target Protein: Human OCT2

**Quantity Size:** 3ml / 6ml (ready to ues)

Clonality: Monoclonal

Clone No.: 5A11

**Isotype:** IgG1, kappa

**Purity:** ≥95%, Purified by Protein G

Form: Liquid

Storage Buffer: 0.01M PBS(pH7.4) with 1% BSA, 0.1% Proclin300.

Storage: Store at 2-8°C for one year. Avoid repeated freeze/thaw cycles.

**Application:** IHC, IF (ready to ues)

**Background:** Oct2 is a transcription factor that specifically binds to the octamer motif (5'-ATTTGCAT-3'). Oct2 regulates transcription in a number of tissues in addition to activating immunoglobulin gene expression. It also modulates transcription transactivation by NR3C1, AR and PGR. Oct2 is B cell specific, belongs to the POU transcription factor family class 2 subfamily and contains 1 homeobox domain.