

**\** 400-901-9800

support@bioss.com.cn

## **Product Datasheet**

Name: Mouse Anti-D-DIMER Monoclonal Antibody

**Description:** Hybridoma clones have been derived from hybridization of myeloma cells with spleed cells of BALB/c mouse immunized with native D-dimer protein derived from human plasma.

Catalog No.	Isotype	Unit	Usage	Buffer
V1102	lgG1	mg	Detection	10mM PBS (pH7.4)

**Specificity:** Mabs react with antigen human D-DIMER

Host: Mouse

Clonality: Monoclonal

Format: Liquid

**Concentration**: ≥1mg/ml

**Purification:** ≥90% (SDS-PAGE)

Preservative: 0.1%Proclin300

**Application:** Recommended for sandwich immunoassays in ELISA and CLIA. Each laboratory

should determine an optimum working titer for use in its particular application.

Best pairs for human D-DIMER complex immunodetection (capture-detection):

Capture antibody: V1101

Detection antibody: V1102

**Storage:** Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Background:** D-DIMER is A fibrin degradation fragment or product that is produced by the action of plasmin on fibrin in the clot dissolution process. See Disseminated intravascular coagulation, Fibrin split products.

**Note:** This product as supplied is intended for research or further manufacturing use only.