

\$400-901-9800
 sales@bioss.com.cn
 techsupport@bioss.com.cn

Rabbit anti-digoxin, whole serum

Cat. Number:	bs-0356Rs
Quantity size:	1ml
Background:	Digoxin is a beta blocker drug originally derived from the foxglove plant, Digitalis lanata. Digoxin is used primarily to improve the pumping ability of the heart in congestive heart failure (CHF), and treat problems such as high blood pressure. It is also used to help normalize some dysrhythmias (abnormal types of heartbeat).
Specificity:	 Rabbit Polyclonal IgG, whole antiserum, unpurified. Reacts with Digoxin Immunogen: KLH conjugated Digoxin
Application:	 Western blotting: 1:10000-50000 Immunohistochemistry: 1:2000-5000 Immunocytochemistry: 1:2000-5000 ELISA: 1: 10000-50000 Optimal working dilutions must be determined by the end user.

Storage: Shipped at 4°C, Store at -20°C (Avoid repeated freeze/thaw cycles).

Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.