

**400-901-9800** 

support@bioss.com.cn

## Rabbit IgG, Isotype Control

Cat. Number: bs-30295P

Quantity Size: 100ul, 500ul

**Purity:** ≥95%, Purified by Protein G

Form: Liquid

**Concentration:** 1ug/ul

**Buffer:** 0.01M PBS, 0.5%BSA, 0.03%Proclin300.

Storage: Store at  $2-8^{\circ}$ C.

**Application:** Isotype Control

0.125ug/1X10<sup>6</sup> cells

Isotype controls are negative controls that are designed to measure the amount of non-specific background signal caused by the primary antibody. Isotype Control have no specificity for the target protein, yet retain all of the non-specific characteristics of an antibody.