

**** 400-901-9800

≥ sales@bioss.com.cn

Product Datasheet

Name: D-dimer Antigen, Native > 90%

Description: Natural Human D-dimer is purified from Human Plasma.

Catalog No.	Unit	Usage	Buffer
bs-41172P	μд	Antigen	Sodium Chloride Buffer, pH 7.5

Host: Human Plasma

Format: Liquid

Concentration: ≥1mg/ml

Purification: ≥90% (SDS-PAGE)

Application: Immunogen, Controls; Recommended for sandwich immunoassays in ELISA and CLIA. Each laboratory should determine an optimum working titer for use in its particular

application.

Storage: Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Background: Cardiac Markers. D-DIMER is A fibrin degradation fragment or product that is produced by the action of plasmin on fibrin in the clot dissolution process.

Note: This product as supplied is intended for research or further manufacturing use only.