

**\** 400-901-9800

support@bioss.com.cn

## **Product Datasheet**

Name: Mouse Anti-SARS-CoV-2 N protein Monoclonal Antibody

**Description:** Hybridoma clones have been derived from hybridization of myeloma cells with

spleen cells of BALB/c mouse immunized with recombinant SARS-CoV-2 N protein.

| Catalog No. | Isotype | Clone No. | Usage              | Buffer           |
|-------------|---------|-----------|--------------------|------------------|
| bsm-41412M  | lgG2b   | 7A7       | Capture /Detection | 10mM PBS (pH7.4) |

**Specificity:** Mab react with recombinant antigen SARS-CoV-2 N protein

Host: Mouse

Clonality: Monoclonal

Format: Liquid

**Concentration:** ≥1 mg/ml

**Purification:** ≥90% (SDS-PAGE)

**Preservative:** 0.1% Proclin300

**Application:** Recommended for sandwich immunoassays in ELISA and CLIA. Each laboratory

should determine an optimum working titer for use in its particular application.

**Storage:** Store at -20 °C for three years. Avoid repeated freeze/thaw cycles.

**Background:** Nucleocapsid protein is a most abundant protein of coronavirus. During virion assembly, N protein binds to viral RNA and leads to formation of the helical nucleocapsid. Nucleocapsid protein is a highly immunogenic phosphoprotein also implicated in viral genome replication and in modulating cell signaling pathways. Because of the conservation of N protein sequence and its strong immunogenicity, the N protein of coronavirus is chosen as a diagnostic tool.

**Note:** This product as supplied is intended for research or further manufacturing use only.