

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
Biomarkers	bs-4708R-FITC	AFP	Hu Monkey
	bs-4708R-HRP	AFP	Hu Monkey
	bs-4821R-Bio	AFP	Hu
	bs-4821R-HRP	AFP	Hu
	bsm-33081M-Bio	AFP	Hu
	bsm-1622M-AP	AFP(A4)	Hu
	bsm-1622M-Cy3	AFP(A4)	Hu
	bsm-1622M-Gold	AFP(A4)	Hu
	bsm-1622M-Bio	AFP(A4)	Hu
	bsm-33053M-HRP	AMACR	Hu Mo Rat
	bs-0091R-AF555	CA125	Hu
	bs-0091R-Cy5	CA125	Hu
	bs-0091R-HRP	CA125	Hu
	bs-0060R-AF488	CEA	Hu Rat
	bsm-4792M-HRP	CEA	Hu
	bsm-4793M-FITC	CEA	Hu
	bsm-1624M-AF488	CEA(B5)	Hu
	bsm-1624M-Bio	CEA(B5)	Hu
	bsm-1624M-Gold	CEA(B5)	Hu
	bsm-1623M-Bio	CEA(C3)	Hu
	bs-6532R-FITC	Epithelial Stromal Interaction 1	Hu Mo Rat Dog Pig Hor Rab Guinea Pig
	bs-0051R-Bio	Melan A	Hu
	bs-19295R-FITC	NMES1	Hu Rat
	bsm-1625M-HRP	PSA(Pa3)	Hu
	bs-0091R-PE	CA125	Hu
	bs-4029R-FITC	Carbonate dehydratase IX	Hu Mo Rat Dog Pig Cow Hor
	bs-0527R-PE	CD10	Hu Mo Rat
	bs-1523R-Bio	CD63	Hu
	bs-1523R-PE-Cy5	CD63	Hu
	bs-1226R-FITC	GPA33	Hu Mo Rat
bs-1226R-PE	GPA33	Hu Mo Rat	
bs-1901R-Bio	STEAP1	Hu Mo Rat Cow Hor	
bs-12850R-PE	SUSD2	Hu Mo Rat	
Extracellular Proteins	bs-5152R-FITC	Alpha Dystroglycan	Hu Mo Rat Chi Dog Pig Hor Rab
	bs-0292R-Cy3	Bovine Serum Albumin	Cow
	bs-0292R-HRP	Bovine Serum Albumin	Cow
	bsm-0945M-PE	Serum albumin (3F4)	Hu Mo Pig
	bs-9005R-FITC	CECR1	Hu Pig
	bs-9005R-PE	CECR1	Hu Pig

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
Extracellular Proteins	bs-1588R-AF555	CK20	Hu Mo Rat
	bs-1588R-FITC	CK20	Hu Mo Rat
	bs-1588R-PE	CK20	Hu Mo Rat
	bs-0709R-AF488	Collagen II	Hu Mo Rat Chi Dog Pig Cow Rab Guinea Pig
	bs-4595R-Cy3	Collagen IV	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-10316R-FITC	DSPP	Hu
	bs-1756R-AF555	Elastin	Hu Mo Rat
	bs-1756R-Bio	Elastin	Hu Mo Rat
	bs-13455R-FITC	Fibronectin/Ugl-Y3	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-1240R-FITC	Hu Fibrinogen	Hu Mo Rat
	bs-1240R-HRP	Hu Fibrinogen	Hu Mo Rat
	bs-15410R-HRP	HAPLN1	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-6894R-HRP	HFREP1	Hu Mo Rat Cow Hor Sheep
	bs-5072R-FITC	HSPG2/Heparan Sulfate Proteoglycan 2	Hu Mo Rat Pig Cow Hor
	bs-7018R-AF555	ITIH1	Hu Mo Rat
	bs-2925R-AF647	Klotho	Hu Mo Dog Pig Cow Hor Rab
	bs-2925R-PE	Klotho	Hu Mo Rat Dog Pig Cow Hor Rab
	bs-18641R-HRP	MAMDC2	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-18693R-FITC	Matrilin 2	Hu Mo Rat Dog Pig Hor
	bs-18694R-FITC	Matrilin 3	Hu Mo Rat Dog Pig Cow Hor Rab Guinea Pig Cat
	bs-7880R-FITC	MEGF5	Hu Mo Rat Chi Pig Cow Rab Sheep
	bs-1704R-Cy5	PAI1	Hu Mo Rat Dog Pig Cow Rab Sheep
	bs-10036R-Cy5	REG4	Hu Mo Rat Dog Pig Cow
	bs-10378R-Cy7	SP-C	Hu Mo Rat Cow Rab Sheep
	bs-0535R-PE	TFF3	Hu Mo Rat Cow
	bs-9622R-FITC	Tropomyosin	Hu Mo Rat Dog Zebrafish Sheep
	bs-2533R-FITC	Versican	Hu Mo Rat Dog Pig Cow Hor
	bs-2533R-RBITC	Versican	Hu Mo Rat Dog Pig Cow Hor
	bs-3574R-Cy3	ADAM10/MADM	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-5856R-HRP	ADAMTS13	Hu Mo Rat Dog Pig Cow Hor

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
Extracellular Proteins	bs-5853R-AF647	ADAM23	Hu Mo Rat Pig Cow Hor Rab
	bs-4605R-FITC	MMP2	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-0423R-HRP	MMP7	Hu Mo Rat
	bs-0397R-AF555	MMP9	Hu Mo Rat
	bs-0397R-FITC	MMP9	Hu Mo Rat
	bs-0397R-RBITC	MMP9	Hu Mo Rat
	bs-3576R-FITC	DTSF/HB EGF	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-3576R-HRP	DTSF/HB EGF	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-2008R-FITC	EGF	Hu
	bs-9762R-FITC	FGF18	Hu Mo Rat Chi Cow Rab Sheep
	bs-9762R-HRP	FGF18	Hu Mo Rat Chi Cow Rab Sheep
	bs-10062R-PE	HGFA Inhibitor 2	Hu Mo Rat Dog Pig Cow Hor Guinea Pig
	bs-0196R-FITC	PDGF-A	Hu Mo Rat Rab
	bs-4938R-Bio	SDF1	Hu Mo Rat Chi Dog Pig Guinea Pig
bs-14001R-HRP	ERPINA6/Cortisol Binding Globulin	Hu Mo Cow Sheep	
Membrane Proteins	bs-5758R-AF488	FAPA	Hu Mo Rat Dog Pig Cow Hor
	bs-16883R-Bio	Laeverin/Aminopeptidase Q	Hu Mo Rat Dog
	bs-16883R-FITC	Laeverin/Aminopeptidase Q	Hu Mo Rat Dog
	bs-0398R-FITC	Annexin V	Hu Mo Rat Dog Pig Cow Guinea Pig G
	bs-0398R-PE	Annexin V	Hu Mo Rat Dog Pig Cow Guinea Pig G
	bs-7669R-FITC	AX2R	Hu Mo
	bs-12418R-Cy3	BAMBI	Hu Mo Rat Dog Pig Cow Rab Sheep Guinea Pig
	bs-8434R-Bio	BTNL9	Hu Mo Rat
	bs-3610R-FITC	Caveolin-3	Hu Mo Rat Dog Cow Hor Rab
	bs-0684R-FITC	CD147	Hu Mo Rat Chi Dog Pig Cow Rab Guinea Pig
	bs-1523R-AF488	CD63	Hu
	bs-1523R-APC	CD64	Hu
bs-1523R-PE	CD65	Hu	

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
Membrane Proteins	bs-4913R-FITC	ENPP1	Hu Mo Rat Dog Pig Hor Rab
	bs-1568R-FITC	ENPP3	Hu Mo Rat Dog Pig Cow
	bs-1568R-PE	ENPP3	Hu Mo Rat Dog Pig Cow
	bs-13372R-Bio	GLEPP1	Hu Mo Rat Dog Hor
	bs-2426R-FITC	Glypican 1	Hu Mo Rat Dog Cow Hor
	bs-13450R-Bio	Glypican 2	Hu Mo Rat
	bs-1112R-AF488	glypican 3	Hu Mo Rat
	bs-1112R-Bio	glypican 3	Hu Mo Rat
	bs-1112R-FITC	glypican 3	Hu Mo Rat
	bs-2177R-FITC	Glypican 6	Hu Mo Rat Dog Hor Rab
	bs-9201R-FITC	IGSF11	Hu Mo Rat Cow Hor Rab Sheep
	bs-16601R-FITC	ILDR2	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-16627R-HRP	IMPG2	Hu Mo Rat Cow Rab Sheep
	bs-18240R-FITC	LHFPL3	Hu Mo Dog Cow Hor Sheep
	bs-5925R-FITC	MAS1	Hu Mo Rat Dog Pig Cow Hor Rab Sheep Guinea Pig
	bs-5925R-PE	MAS1	Hu Mo Rat Dog Pig Cow Hor Rab Sheep Guinea Pig
	bs-6454R-FITC	MFGE8/BA46	Hu Rab
	bs-10673R-AF647	Podoplanin	Hu Mo Rat
	bs-20015R-APC	Placental alkaline phosphatase (PLAP)	Hu Mo
	bs-17482R-Bio	Sialidase 3	Hu Dog Cow
	bs-1770R-FITC	Spastin	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-1309R-PE	Syndecan-1	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-2101R-Cy5	TEM1	Hu Mo Rat Dog Pig Cow Hor
	bs-2101R-FITC	TEM1	Hu Mo Rat Dog Pig Cow Hor
	bs-2101R-PE	TEM1	Hu Mo Rat Dog Pig Cow Hor
	bs-6241R-FITC	TM4SF3	Hu Mo Rat

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
	bs-9410R-FITC	TSPAN1	Hu
	bs-10025R-RBITC	Ulex Europaeus Lectin 1	Ulex europaeus
	bs-3584R-Bio	Ulex Europaeus Lectin 1	Ulexeuropaeus
	bs-18479R-HRP	ZDHHC19	Hu
Channels	bs-1604R-FITC	ABCB5	Hu Mo Rat Dog Cow Hor
	bs-1253R-FITC	AQP3	Hu Mo Rat
	bs-4611R-Cy5	Aquaporine 2	Hu Mo Rat Dog
	bs-4611R-PE	Aquaporine 2	Hu Mo Rat Dog
Channels	bs-20129R-FITC	SCN1A	Hu Mo Rat Dog Cow Hor Rab Sheep
Receptors	bs-4119R-Cy3	ASGPR1	Hu Mo Rat Pig Cow Hor Rab
	bs-4119R-FITC	ASGPR1	Hu Mo Rat Pig Cow Hor Rab
	bs-4119R-HRP	ASGPR1	Hu Mo Rat Pig Cow Hor Rab
	bs-4119R-PE	ASGPR1	Hu Mo Rat Pig Cow Hor Rab
	bs-4825R-Bio	ASGPR1	Hu Mo Rat Hor
	bs-1860R-FITC	CALCRL	Hu Mo Rat
	bs-4194R-FITC	CD167b/DDR2	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-4736R-FITC	CD168	Hu Mo Rat Cow Hor
	bs-0988R-AF488	CD71/Transferrin receptor	Hu Mo Rat Dog Pig Cow Hor Guinea Pig
	bs-1782R-APC	CD71/Transferrin receptor	Hu Mo Rat Chi Dog Pig Cow Hor
	bs-1782R-Bio	CD71/Transferrin receptor	Hu Mo Rat Chi Dog Pig Cow Hor
	bs-2052R-RBITC	Transferrin	Hu
	bs-20716R-PE	c-Kit/CD117	Hu Mo Hor
	bs-20717R-APC	c-Kit/CD117	Hu Mo Rat Dog Hor
	bs-14284R-HRP	Dexas1	Hu Mo Rat Chi Dog Pig Cow Hor Rab Sheep
	bsm-33140M-Cy5	EGFR	Hu
	bsm-33140M-FITC	EGFR	Hu
	bsm-33050M-FITC	EGFR	Hu Mo Rat Dog Cow Sheep
	bs-0165R-PE	EGFR	Hu Mo Rat Dog Pig
	bs-10007R-AP	EGFR	Hu Mo Rat Pig Cow Rab Sheep
bs-2558R-AF647	EGFRvIII	Hu	

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
Receptors	bs-10627R-FITC	EGFRvIII	Hu Mo Rat Dog Pig Hor Rab Sheep
	bs-1454R-FITC	ErbB 3	Hu Mo Rat Chi Dog Cow Hor Rab
	bsm-2156M-HRP	HER2 receptor	Hu Dog Hor
	bs-3021R-FITC	phospho-ALK (Tyr1278 + Tyr1282 + Tyr1283)	Hu Mo Rat Chi Dog Cow Hor
	bs-13170R-HRP	phospho-FGFR3 (Tyr724)	Hu Mo Rat Dog Cow Hor Rab Sheep
	bs-2674R-FITC	phospho-VEGFR2 (Tyr951)	Hu Mo Rat Rab
	bs-2202R-FITC	VEGFR3	Hu Mo Rat
	bs-2202R-PE	VEGFR3	Hu Mo Rat
	bs-13508R-FITC	GPCR G2A	Hu
	bs-13513R-PE	GPR110	Hu Mo
	bs-16292R-FITC	GPR156	Hu Mo Rat
	bs-15394R-FITC	GPRC5C	Hu Mo Rat Pig Sheep
	bs-10368R-HRP	LPA2	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-10368R-PE	LPA2	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-1419R-FITC	MC1 Receptor	Hu Mo Rat
	bs-1419R-PE	MC1 Receptor	Hu Mo Rat
	bs-12699R-FITC	PLA2R	Hu Mo Cow Rab
	bs-4690R-FITC	Tissue factor/CD142	Hu Mo Rat Dog Pig Cow Hor Rab Guinea Pig
bs-4690R-PE	Tissue factor/CD142	Hu Mo Rat Dog Pig Cow Hor Rab Guinea Pig	
	bs-2016R-AF488	Integrin Alpha V + Beta 1	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-2016R-Bio	Integrin Alpha V + Beta 1	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-2016R-Cy3	Integrin Alpha V + Beta 1	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-2016R-HRP	Integrin Alpha V + Beta 1	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-2016R-PE	Integrin Alpha V + Beta 1	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-2016R-PE-Cy7	Integrin Alpha V + Beta 1	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-1356R-Cy5	Integrin Alpha V + Beta 5	Hu Mo Rat
	bs-2095R-PE	Integrin alpha 1	Hu Mo Rat Cow Hor

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
Adehesion and Junction	bs-0641R-Cy7	Integrin alpha 4	Hu Mo Rat Cow Hor Rab
	bs-2055R-Bio	Integrin alpha 4+Beta 7	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-0342R-AF594	Integrin beta 3	Hu Mo Pig Cow Rab
	bs-0342R-AF647	Integrin beta 3	Hu Mo Pig Cow Rab
	bs-16615R-PE-Cy5	Integrin beta 8	Hu Mo Rat Chi Pig Rab
	bs-13771R-FITC	ITGA11	Hu Mo Rat Chi Dog Cow Hor Rab
	bs-0901R-PE	LAMR1	Hu Mo Rat Dog
	bs-0300R-FITC	Mesothelin	Hu
	bs-20623R-Bio	N-cadherin	Hu Mo Rat Pig Cow Hor Rab Sheep
	bs-1159R-FITC	P-cadherin	Hu Mo Rat Chi Dog Cow Rab Guinea Pig
	bs-11109R-FITC	PCDH12	Hu Mo Rat Cow Hor Sheep
	bs-11110R-FITC	PCDH17	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-11110R-HRP	PCDH17	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-20154R-FITC	PCDH7	Hu Mo Rat Pig Cow Sheep
	bs-11114R-FITC	PCDH8	Hu Mo Rat Pig Cow Rab Goose Sheep Chimpanzee
	bs-11115R-FITC	PCDHA1	Hu Mo Rat Dog
	bs-6026R-FITC	TSLC1	Hu Mo Rat Dog Pig Cow Hor Rab
	bs-6026R-PE	TSLC1	Hu Mo Rat Dog Pig Cow Hor Rab
	bs-0396R-Cy5.5	VCAM-1	Hu Mo Rat
	bs-0878R-Cy3	CD144/VE Cadherin	Hu Mo Rat
bs-0878R-FITC	CD144/VE Cadherin	Hu Mo Rat	
bs-0878R-PE	CD144/VE Cadherin	Hu Mo Rat	
Adehesion and Junction	bs-1618R-Bio	CD146	Hu Mo Rat Dog Cow Hor
	bs-1618R-PE	CD146	Hu Mo Rat Dog Cow Hor
	bs-4736R-PE	CD168	Hu Mo Rat Cow Hor
	bs-0503R-FITC	CD18	Hu Mo Rat Pig Cow Hor Sheep

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
	bs-0593R-Bio	EpCAM	Hu Mo Rat Dog Pig Cow Rab
	bs-1513R-AF647	EpCAM	Hu Mo Rat
	bs-1513R-Bio	EpCAM	Hu Mo Rat
	bs-1513R-Cy7	EpCAM	Hu Mo Rat
	bs-1513R-FITC	EpCAM	Hu Mo Rat
	bs-1513R-PE	EpCAM	Hu Mo Rat
	bs-16158R-Bio	Focadhesin	Hu Mo Rat Dog Pig
	bs-2383R-AF647	GJB3	Hu Mo Rat Dog Pig Cow Hor Rab Guinea Pig
	bs-2788R-FITC	Claudin 3	Hu Mo Rat Dog Pig Cow Sheep
bs-2790R-Bio	Claudin 4	Hu Mo Rat Chi Dog Pig Cow Hor Rab Shee	
Transporters	bs-0657R-FITC	MRP1/ABCC1	Hu Mo Rat Chi Cow
	bs-10357R-FITC	SLC22A12	Hu
	bs-1958R-APC	SLC10A1	u Mo Rat Dog Hor Rab
	bs-1958R-PE	SLC10A1	u Mo Rat Dog Hor Rab
	bs-19804R-FITC	SLC20A1	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-21248R-HRP	SLC26A9	Hu Mo Rat
	bs-21278R-FITC	SLC38A5	Hu Mo Rat
bs-19826R-FITC	SLC38A10	Hu Mo Rat Dog Hor	
Vesicle Trafficking	bs-2932R-AF647	Clathrin heavy chain	Hu Mo Rat Dog Pig Cow Rab Bee bee ant frog silkworm
	bs-2932R-Bio	Clathrin heavy chain	Hu Mo Rat Dog Pig Cow Rab Bee bee ant frog silkworm
Cytoskeletal Proteins	bsm-33039M-FITC	alpha Tubulin	Hu Mo Rat Chi Dog Pig Rab Guinea Pig Danio rerio
	bs-0159R-FITC	Tubulin-alpha	Hu Mo Rat Cow Sheep
	bsm-33041M-AF555	beta I Tubulin	Hu Mo Rat
	bsm-33041M-PE-Cy5	beta I Tubulin	Hu Mo Rat
	bs-1060R-PE	Cytokeratin 5	Hu Mo Rat Dog Pig Cow Hor Rab
	bs-1106R-FITC	Cytokeratin 8	Hu Mo Rat Chi Dog Pig
	bs-4625R-AF488	CK12	Hu Mo Rat Dog Rab
	bs-4625R-FITC	CK12	Hu Mo Rat Dog Rab
	bsm-33062M-PE	CK16	Hu Mo Rat
	bs-1028R-AF488	Cytokeratin 19	Hu Rat
bs-1026R-FITC	Desmin	Hu Mo Rat Chi Dog Pig Cow Hor Rab Sheep Guinea Pig	

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
Cytoskeletal Proteins	bs-14360R-FITC	DNAH11	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-14363R-Bio	DNAH2	Hu Mo Rat Pig Hor
	bs-11024R-FITC	DNAI1 + 2	Hu Mo Rat Dog Cow Hor Rab Sheep
	bs-5885R-AF555	Fast Myosin Skeletal Heavy chain	Hu Mo Rat Chi Dog Pig
	bs-1160R-Bio	Gelsolin	Hu Mo Rat
	bs-11028R-FITC	MAP1B	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-1534R-FITC	MAP1LC3A	Hu Mo Rat Chi Pig Cow
	bs-11513R-FITC	MKS1	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-6075R-FITC	MTSS1	Hu Mo Rat Chi Dog Cow Hor Rab
	bs-9862R-Cy7	MYH7B	Hu Mo Rat Chi Hor
	bs-10366R-FITC	MYPN	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-19209R-FITC	Nesprin 3	Hu Mo Rat
	bs-3539R-FITC	Paxillin	Hu Mo Rat Dog Cow Hor Rab
	bs-10367R-Cy3	Sarcomeric Actinin Alpha	Hu Mo Rat Chi Dog Pig Cow Hor Rab Sheep
Signal Transduction	bs-3916R-FITC	ADCY10	Hu Mo Rat Dog Pig Cow Hor Rab
	bs-1332R-FITC	Beta arrestin 2	Hu Mo Rat Chi Dog Pig Cow Hor Rab Guinea Pig
	bs-14404R-FITC	DOCK8	Hu Mo Rat Dog Cow Hor Sheep
	bs-16017R-FITC	FAM83G	Hu Mo Rat Cow Hor Rab
	bs-13302R-PE-Cy7	GBP1	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-13303R-FITC	GBP4	Hu Mo Rat Pig Cow Hor Rab Sheep
	bs-13249R-Bio	GNAI1	Hu Mo Rat Pig Cow Zeb Sheep
	bs-7844R-FITC	GSK3B interacting protein/Ninein	Hu Mo Rat Dog Pig Cow Hor Rab
	bs-19108R-FITC	MURC	Hu Mo Rat Hor Rab Sheep
	bs-5590R-FITC	phospho-14-3-3 epsilon (Thr232)	Hu Mo Rat Cow Sheep

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
Signal Transduction	bs-3341R-FITC	Phospho-PLC beta3 (Ser1105)	Hu Mo Rat Dog Pig Cow Hor
	bs-18843R-PE	RRAGD	Hu Mo Rat Pig Cow Hor Rab Sheep
	bs-12941R-HRP	CSK	Hu Mo Rat Chi Dog Pig Cow Sheep Monkey
	bs-1687R-FITC	DAPK1	Hu Mo Rat
	bs-2768R-PE	EIF2AK4/GCN2	Hu Mo Rat Dog Hor Rab
	bs-14458R-HRP	Inactive DUSP27	Hu Mo Rat Cow Hor Rab Sheep
	bs-4134R-FITC	MAP4K1	Hu Mo Rat Cow
	bs-5224R-HRP	phospho-B-Raf (Thr401)	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-3155R-FITC	Phospho-GCN2 (Thr899)	Hu Mo Rat Dog Cow Hor
	bs-3155R-PE	Phospho-GCN2 (Thr899)	Hu Mo Rat Dog Cow Hor
	bs-5491R-FITC	phospho-MAP4K4 (Ser629)	Hu Mo Rat Chi Dog Hor Rab
	bs-5492R-FITC	phospho-MAP4K4 (Ser631)	Hu Mo Rat Chi Dog Hor Rab
	bs-18672R-FITC	Phospho-MAPK6 (Ser189)	Hu Mo Rat Chi Dog Pig Cow Hor Rab Sheep
	bs-3280R-FITC	Phospho-MLK3 (Thr277 + Ser281)	Hu Mo Rat Dog Pig Cow Rab Guinea Pig
	bs-3280R-HRP	Phospho-MLK3 (Thr277 + Ser281)	Hu Mo Rat Dog Pig Cow Rab Guinea Pig
	bs-5570R-HRP	phospho-PI3KCA (Tyr317)	Hu Mo Rat
	bs-11327R-FITC	PTP zeta	Hu Mo Rat
	bs-11327R-PE	PTP zeta	Hu Mo Rat
	bs-4904R-PE	SGK1	Hu Mo Rat Dog Cow Hor
	bs-0059R-FITC	APAF1(NT)	Hu Mo Rat Chi Dog Cow Hor
	bs-0738R-HRP	Livin	Hu Mo Rat Dog Cow
	bsm-33339M-FITC	MLKL	Hu
bs-0037M-FITC	AIF	Hu Mo Rat Chi Dog Pig Cow Rab Sheep	
bs-0052R-FITC	Caspase 8	Hu Mo Rat Chi Dog Pig Cow Sheep	
bs-2322R-AF488	Caspase-7 p11	Hu Mo Rat Chi Dog Pig Cow Hor Rab	

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
Signal Transduction	bs-1281R-FITC	XIAP/BIRC4	Hu Mo Rat
	bs-10442R-APC	Caspase-1 p20	Hu Mo Rat
	bs-10442R-FITC	Caspase-1 p20	Hu Mo Rat
	bs-9074R-APC	CDC37L1	Hu Mo Rat Dog Pig Cow Hor Sheep
Cytosolic Proteins	bs-1611R-Bio	Cathepsin K	Hu Mo Rat Dog Pig Cow Hor Guinea Pig
	bs-1611R-PE-Cy7	Cathepsin K	Hu Mo Rat Dog Pig Cow Hor Guinea Pig
	bs-13888R-HRP	CHCHD10	Hu Mo Rat Pig Cow Sheep
	bs-4100R-FITC	DRP1	Hu Mo Rat Dog Hor
	bs-14991R-FITC	FAM38A	Hu Mo Rat Dog Cow Sheep
	bs-16010R-FITC	FAM73A	Hu Mo Rat
	bs-8155R-AF488	GM130	Hu Mo Rat Pig Cow
	bs-2556R-FITC	LAMP3	Hu Mo Rat Cow Hor Rab
	bs-1970R-AF488	LAMPL1	Hu Mo Rat
	bs-1970R-Bio	LAMPL1	Hu Mo Rat
	bs-1970R-FITC	LAMPL1	Hu Mo Rat
	bs-7357R-AF647	TOMM20	Hu Mo Rat Dog Cow Hor Rab Sheep
	bs-13922R-FITC	Chorein	Hu Mo Rat Cow Hor Sheep
	bs-18862R-AF647	RRS1	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-8034R-AF488	COMMD1	Hu Mo Rat Cow Hor Rab Sheep
	bs-2787R-FITC	LMP2	Hu Mo Rat Dog Pig Cow Hor Rab
	bs-9353R-FITC	Proteasome 20S alpha 3(Ser250)	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-13106R-Cy3	ERp72	Hu Mo Rat Cow Rab Sheep Monkey
	bs-5913R-AF488	Calreticulin	Hu Mo Rat Chi Dog Pig Cow Hor Rab Guinea Pig
	bs-5913R-FITC	Calreticulin	Hu Mo Rat Chi Dog Pig Cow Hor Rab Guinea Pig
	bs-5913R-PE	Calreticulin	Hu Mo Rat Chi Dog Pig Cow Hor Rab Guinea Pig
	bs-17865R-FITC	MT-ATP6	Hu
bs-3888R-Bio	NDUFS3	Hu Mo Rat Chi Dog Pig Cow Hor Rab	

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
Cytosolic Proteins	bs-9601R-Bio	Frataxin	Hu Mo Rat Chi Dog Pig Cow Hor
	bs-19867R-HRP	RHBDF1	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-16870R-FITC	KMO	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-16870R-PE	KMO	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-1061R-AF594	Myeloperoxidase	Hu Mo Rat Dog Cow Hor Rab Guinea Pig
	bs-1061R-FITC	Myeloperoxidase	Hu Mo Rat Dog Cow Hor Rab Guinea Pig
	bs-24004R-Bio	NEU2	Mo Rat
	bs-12788R-PE-Cy3	TRAP 5	Hu Mo Rat
	bs-14019R-FITC	CPEB2	Hu Mo Rat Pig Cow Hor Sheep
	bs-16709R-FITC	TRBP	Hu Mo Rat Pig Cow Hor Sheep
	bs-16375R-FITC	GW182	Hu Mo Rat
	bs-17346R-AF555	hnRNPLL	Hu Mo Rat Chi Dog Pig Cow Hor Rab Sheep
	bs-2761R-APC	DBC1	Hu Mo Rat Chi Pig Cow Hor Rab
	bs-16015R-HRP	FAM83B	Hu Mo Rat Rab
	bs-16018R-FITC	FAM83H	Hu Mo Rat Dog Cow Hor Sheep
	bs-3657R-FITC	FLASH	Hu Mo Rat Dog Pig Cow Hor
	bs-23062R-PE	Involucrin	Hu Mo Rat
	bs-8065R-FITC	PID1	Hu Mo Chi Rab
bs-3578R-HRP	TCTP	Hu Mo Rat Dog Pig Cow Hor Rab Sheep	
Nuclear Proteins, Transcription Regulation	bs-9092R-AF594	hnRNP U	Hu Mo Rat Pig Cow
	bs-3477R-AF594	Phospho-YB1 (Ser102)	Hu Mo Rat Chi Dog Pig Cow Hor Rab Guin
	bs-19512R-FITC	NUDT16	Hu Rat
	bs-19597R-AF555	PUS7	Hu Mo Rat Chi Pig Rab Sheep
	bs-17000R-HRP	G2E3	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-2935R-FITC	Nucleoporin p62	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-19340R-FITC	NR2F6	Hu Mo Rat Cow Hor Sheep
	bs-1485R-HRP	BIN1	Hu Mo Rat Dog Cow Hor
	bs-16676R-FITC	INSM2	Hu Mo Rat Pig Cow Rab Sheep

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
	bs-17090R-HRP	KBTBD7	Hu Mo Rat Dog Pig Cow Hor Rab Sheep
	bs-23199R-FITC	MALAT1	Hu
	bs-23199R-PE	MALAT1	Hu
	bs-8536R-FITC	Nucleolin/C23	Hu Mo Rat Dog Pig Cow Hor
	bs-2090R-AF647	P53(wt-p53)	Hu Mo Rat Dog Pig Hor Rab Sheep
	bs-0723R-FITC	P63	Hu Mo Rat Dog Pig Cow Hor Rab Sheep Guinea Pig
	bs-3582R-FITC	TIF1 gamma	Hu Mo Rat Chi Dog Pig Cow Hor
	bs-12809R-FITC	ZNF385B	Hu Mo Rat Dog Pig Cow Sheep
	bs-8622R-FITC	KPNA2	Hu Mo Rat Pig Cow Hor Rab
Cell Cycle, DNA Damage	bs-0633R-FITC	CDK4	Hu Mo Rat Pig Cow
	bs-0740R-AF647	CDKN2A/p16-INK4a	Hu Mo Rat Pig
	bs-0740R-FITC	CDKN2A/p16-INK4a	Hu Mo Rat Pig
	bs-1174R-Cy5	CDKN2A/p19ARF	Hu Mo Rat
	bs-1174R-FITC	CDKN2A/p19ARF	Hu Mo Rat
	bsm-33321M-FITC	Cyclin B1	Hu Mo Rat
	bs-0660R-FITC	Cyclin D3	Hu Mo Rat Dog Pig Cow Rab
	bs-18952R-FITC	MLIP/C6orf142	Hu Mo Rat Dog Cow Sheep
	bs-17851R-FITC	MSH4	Hu Mo Rat Pig Cow Hor Sheep
	bs-7711R-FITC	PCNT	Hu Mo Rat Chi Dog Pig Cow Hor Rab
	bs-3293R-FITC	Phospho-Mre11 (Ser676)	Hu Mo Rat Dog
	bs-9173R-APC	RNF8	Hu Mo Rat Dog Pig Hor
	bs-11233R-PE	RPA70	Hu Mo Rat Dog Pig Rab
	bs-17269R-FITC	SCAPER	Hu Mo Rat Chi Dog Pig Cow Hor Rab Sheep
	bs-7360R-FITC	Topoisomerase III	Hu Mo Rat
Cell Division	bs-1962R-PE	ACD	Hu Mo Rat Dog Cow
	bs-0489R-HRP	BrdU	Rat BrdU
	bs-15414R-Cy3	HAUS3	Hu Mo Rat
	bs-23102R-Cy5	Ki-67	Hu Mo
	bs-23102R-FITC	Ki-67	Hu Mo
	bsm-33070M-FITC	Ki-67	Hu Mo Rat
	bs-3512R-RBITC	NuMA	Hu Mo Rat Dog

Cell Biology

(请核对选择抗体所带有的标记物是否适用于您将开展的实验: 1. HRP/Biotin 或 AP 标记的抗体主要推荐用于 WB; 2. 荧光素标记的抗体可用于 IF 或 FCM)

Category	Cat.No.	Target Name	Species of reactivity
	bs-3066R-Cy3	Phospho-CENPA (Ser7)	Hu Rab
	bs-19361R-HRP	SAC3D1	Hu Mo Rat
Protein Modification	bs-8347R-HRP	UBE2I/UBC9	Hu Mo Rat
	bs-7338R-HRP	Sumo 2+3	Hu Mo Rat Chi Cow Hor Rab Zeb
	bs-7621R-FITC	ART4	Hu Mo Rat Pig Cow Rab Sheep
Others	bs-2231R-FITC	LYRIC	Hu Mo Rat Dog Pig Cow
	bs-14726R-HRP	FAM135A	Hu Mo Rat Dog Pig
	bs-14737R-Bio	FAM160A1	Hu Mo Rat Dog Pig Cow Hor Sheep
	bs-14764R-HRP	FAM177A1	Hu Mo Rat Dog Pig
	bs-16983R-Bio	KIAA0802	Hu Mo Rat Chi Dog
	bs-2039R-FITC	GRGDS	Hu Gly-Arg-Gly-Asp-Ser