

ELISA 实验技术指南



公司介绍

博奥森：专注科研，以优质抗体引领未来

博奥森生物，自 2001 年诞生以来，始终坚守在生命科学前沿，为全球科研人员提供卓越品质的免疫学试剂产品与服务。我们拥有资深的科学家团队、先进的抗体发现、验证与生产平台，始终坚持“自主研发、原始创新”的理念，确保每一款产品都能达到国际标准，持续提供“4R”品质科研工具【Repeatable（可重复）、Replicable（可复制）、Reproducible（可再现）和 Reliable（可靠的）】。

博奥森承诺：中国制造，国际品质


经过多年深耕，博奥森已拥有成熟的蛋白表达、抗体制备、抗体验证、单 B 细胞筛选、噬菌体展示、重组抗体及工业化生产等平台，为我们的抗体开发奠定了扎实的基础。

立足本土，面向世界，依托我们位于波士顿的 Bioss 全球品牌中心，携手全球渠道伙伴，我们的产品已遍布全球众多高校、科研机构及生物制药公司，成为科研人员的得力助手。我们坚守“中国制造，国际品质”的承诺，并严格执行 ISO 9001, ISO 13485 国际质量管理体系，以行业高标准为全球生命科学研究提供全方位的服务。

在不断督促自身突破创新，品质为先的基础上，我们尤其珍视每一位客户的反馈，用户的声音是我们不断优化产品、提升品质的重要动力。

博奥森，以创新为魂，以品质为根，与您携手共创生命科学的美好未来。

产品体系




目录

01 ELISA 技术	P02
02 基本原理	P03
03 操作流程	P04
04 常见 ELISA 检测类型	P05
05 ELISA 样品制备方法	P09
06 产品目录	P10
07 ELISA 配套试剂	P51
08 ELISA 试剂盒定制	P52
09 ELISA 常见问题解答指南	P53

01 / ELISA 技术

ELISA 是一种灵敏度（敏感性）高，特异性强，重复性好的实验诊断方法。其将抗原与抗体免疫反应的特异性和酶的高效催化作用原理有机地结合起来，可敏感地测定体液中微量的特异性抗体或抗原。此项技术自 70 年代初问世以来，发展十分迅速，由于其试剂稳定、易保存，操作简便，结果判断较客观等因素，目前已被广泛用于生物学和医学科学的许多领域。




ELISA 是酶联免疫吸附剂测定 (Enzyme-Linked ImmunoSorbent Assay) 的简称。



02 / 基本原理


ELISA 的基础是抗原或抗体的固相化及抗原或抗体的酶标记。结合在固相载体表面的抗原或抗体仍保持其免疫学活性，酶标记的抗原或抗体既保留其免疫学活性，又保留酶的活性。在测定时，受检标本（测定其中的抗体或抗原）与固相载体表面的抗原或抗体起反应。用洗涤的方法使固相载体上形成的抗原抗体复合物与液体中的其他物质分开。再加入酶标记的抗原或抗体，也通过反应而结合在固相载体上。此时固相上的酶量与标本中受检物质的量呈一定的比例。加入酶反应的底物后，底物被酶催化成为有色产物，产物的量与标本中受检物质的量直接相关，故可根据呈色的深浅进行定性或定量分析。由于酶的催化效率很高，间接地放大了免疫反应的结果，使测定方法达到很高的敏感度。



03 / 操作流程

ELISA 主要包括四个大致流程：1、将抗原或抗体包被于聚苯乙烯材质的酶标板上；2、使用不与包被的抗原（抗体）或检测抗原（抗体）发生反应的封闭试剂封闭酶标板；3、加入受检抗原（抗体）和酶标记物发生免疫反应；4、加入底物显色。

ELISA 流程示意图如下：



04 / 常见的 ELISA 检测类型

ELISA 可用于测定抗原，也可用于测定抗体。在这种测定方法中有三个必要的试剂：

- 固相的抗原或抗体，即“免疫吸附剂” (Immunosorbent)；
- 酶标记的抗原或抗体，称为“结合物” (Conjugate)；
- 酶反应的底物。


根据试剂的来源和标本的情况以及检测的具体条件，可设计出各种不同类型的检测方法。用于临床检验的 ELISA 主要有以下几种类型：

直接法测抗原

将抗原直接吸附（或通过捕获抗体结合）到固相载体表面，洗涤除去未结合的抗原及杂质。加入酶标抗体进行孵育，使酶标抗体与固相载体上的抗原结合，形成固相免疫复合物。彻底洗涤除去未结合的酶标抗体。加入底物，酶催化底物反应生成可检测的产物（例如显色或发光产物）。产物的颜色或荧光强度与标本中受检抗原的量成正比，通过测定产物信号即可定量抗原。

操作步骤如下：

- Step-1** 将特异性抗原与固相载体联结，形成固相抗原。洗涤除去未结合的抗原及杂质。
- Step-2** 加酶标抗体，保温反应。酶标抗体与固相抗原结合，形成固相抗原抗体复合物。洗涤除去其他未结合物质。
- Step-3** 加底物显色。固相上的酶催化底物成为有色产物。通过比色，测知标本中抗原的量。




间接法测抗体

间接法是检测抗体常用的方法。其原理为利用酶标记的抗抗体（例如：抗人免疫球蛋白抗体）以检测与固相抗原结合的受检抗体，故称为间接法。

操作步骤如下：

- Step-1** 将特异性抗原与固相载体联结，形成固相抗原。洗涤除去未结合的抗原及杂质。
- Step-2** 加稀释的受检血清，保温反应。血清中的特异抗体与固相抗原结合，形成固相抗原抗体复合物。经洗涤后，固相载体上只留下特异性抗体，血清中的其他成份在洗涤过程中被洗去。
- Step-3** 加酶标抗抗体。可用酶标抗人 IgG 以检测总抗体，但一般多用酶标抗人 IgG 检测 IgG 抗体。固相免疫复合物中的抗体与酶标抗体结合，从而间接地标记上酶。洗涤后，固相载体上的酶量与标本中受检抗体的量正相关。
- Step-4** 加底物显色。



本法主要用于对病原体抗体的检测而进行传染病的诊断。

间接法的优点是只要变换包被抗原就可利用同一酶标抗抗体建立检测相应抗体的方法。

间接法成功的关键在于抗原的纯度。虽然有时用粗提抗原包被也能取得实际有效的结果，但应尽可能予以纯化，以提高试验的特异性。


特别注意除去能与一般健康人血清发生反应的杂质，例如以 *E.coli* 为工程酶的重组抗原，如其中含有 *E.coli* 成份，很可能与受 *E.coli* 感染者血清中的抗 *E.coli* 抗体发生反应。

抗原中也不能含有与酶标抗人 IgG 反应的物质，例如来自人血浆或人体组织的抗原，如不将其中的 IgG 去除，试验中也会发生假阳性反应。另外如抗原中含有无关蛋白，也会因竞争吸附而影响包被效果。

双抗体夹心法测抗原

双抗体夹心法是检测抗原最常用的方法，操作步骤如下：

- Step-1** 将特异性抗体与固相载体联结，形成固相抗体。洗涤除去未结合的抗体及杂质。
- Step-2** 加受检标本，保温反应。标本中的抗原与固相抗体结合，形成固相抗原抗体复合物。洗涤除去其他未结合物质。
- Step-3** 加酶标抗体，保温反应。固相免疫复合物上的抗原与酶标抗体结合。彻底洗涤未结合的酶标抗体。此时固相载体上带有的酶量与标本中受检抗原的量相关。
- Step-4** 加底物显色。固相上的酶催化底物成为有色产物。通过比色，测知标本中抗原的量。



在临床检验中，此法适用于检验各种蛋白质等大分子抗原，例如 HBsAg、HBeAg、AFP、HCG 等。

只要获得针对受检抗原的特异性抗体，就可用于包被固相载体和制备酶结合物而建立此法。


如抗体的来源为抗血清，包被和酶标用的抗体最好分别取自不同种属的动物。

如应用单克隆抗体，一般选择两个针对抗原上不同决定簇的单抗，分别用于包被固相载体和制备酶结合物。

这种双位点夹心法具有很高的特异性，而且可以将受检标本和酶标抗体一起保温反应，作一步检测。

竞争法测抗原

小分子抗原或半抗原因缺乏可作夹心法的两个以上的位点，因此不能用双抗体夹心法进行测定，可以采用竞争法模式。其原理是标本中的抗原和一定量的酶标抗原竞争与固相抗体结合。标本中抗原含量愈多，结合在固相上的酶标抗原愈少，最后的显色也愈浅。小分子激素、药物等 ELISA 测定多用此法。



ABS-ELISA 法

ABS-ELISA 法是一种基于亲和素 - 生物素系统 (ABS) 的酶联免疫吸附测定法 (ELISA)，通过利用亲和素与生物素之间的高亲和力来提高检测的灵敏度。ABS-ELISA 法可以分为两种类型：酶标记亲和素 - 生物素 (LAB) 法和桥联亲和素 - 生物素 (ABC) 法。

ABS-ELISA 法的核心在于亲和素和生物素之间的强结合力和高稳定性。亲和素是一种糖蛋白，每个分子由 4 个能和生物素结合的亚基组成，而生物素为小分子化合物。通过化学方法制成的衍生物生物素 - 羟基琥珀酰亚胺酯可以与抗体或酶形成生物素标记产物，这种标记方法简便且不影响抗体或酶的生物学活性。在检测过程中，特异性抗体被生物素化，酶分子标记在亲和素分子上，通过亲和素与生物素的结合，将酶分子联结到抗原分子上，从而进行酶促反应检测抗原。由于 ABS-ELISA 较普通 ELISA 多用了两种试剂，增加了操作步骤，在临床检验中 ABS-ELISA 应用不多。

05 / ELISA 样品制备方法

样品制备是进行 ELISA 之前的一个重要环节，为方便实验者使用，在此给大家提供一个样品制备通用流程，主要过程如下：

【细胞培养上清】

将细胞培养基移至无菌离心管，在 4℃ 条件下 1500 rpm 离心 10 min，然后将上清等量分装于小 EP 管并于 -80℃ 下保存，避免反复冻融。

【血清】

室温血液自然凝固 20 min 后，在 4℃ 条件下 3000 rpm 离心 10 min，然后将上清等量分装于小 EP 管并于 -80℃ 下保存，保存过程中如有沉淀，请再次离心，避免反复冻融。

【血浆】

将全血收集到含抗凝剂的管中，根据标本的实际要求选择 EDTA、柠檬酸钠或肝素作为抗凝剂，混合 10-20 min 后，在 4℃ 条件下 3000 rpm 离心 10 min，然后将上清等量分装于小 EP 管并于 -80℃ 下保存，避免反复冻融。

【尿液】

用无菌管收集尿液，在 4℃ 条件下 12,000 rpm 离心 5 min，然后将上清等量分装于小 EP 管并于 -80℃ 下保存，避免反复冻融。

【唾液】

用无菌管收集唾液，在 4℃ 条件下 12,000 rpm 离心 2 min，然后将上清等量分装于小 EP 管并于 -80℃ 下保存，避免反复冻融。

注：

在实际样品制备时，应查阅相关文献或材料以根据靶标物的实际情况来建立制备方法；在收集样品前，要做好完整的计划，必须充分了解蛋白是否足够稳定，若不稳定，我们建议新鲜样品最好及时检测，若蛋白稳定又不能及时检测，可以储存于 4℃ 备用，但不要超过 7 天；因冰箱与室温存在一定温差，蛋白在反复冻融时容易降解，直接影响实验质量，所以要尽量避免此步骤。

06 / 产品目录

产品名称	产品编号	产品规格
Human		
Human 5-Hydroxytryptamine Receptor 1A (HTR1A) ELISA Kit	BSKH63694	96T
Human 5-Hydroxytryptamine Receptor 2A (HTR2A) ELISA Kit	BSKH63695	96T
Human A Disintegrin And Metalloproteinase With Thrombospondin 4 (ADAMTS4) ELISA Kit	BSKH63052	96T
Human a1AGP ELISA Kit	BSKH60865	96T
Human a1M ELISA Kit	BSKH60866	96T
Human a2M ELISA Kit	BSKH60371	96T
Human a2ML1 ELISA Kit	BSKH60867	96T
Human a2PI ELISA Kit	BSKH60868	96T
Human Ab1-40 ELISA Kit	BSKH60869	96T
Human Ab1-42 ELISA Kit	BSKH60870	96T
Human ABCG2 ELISA Kit	BSKH60871	96T
Human ACE ELISA Kit	BSKH60372	96T
Human ACE2 ELISA Kit	BSKH60121	96T
Human ACP5 ELISA Kit	BSKH60872	96T
Human ACP5 ELISA Kit	BSKH61337	96T
Human ACTH ELISA Kit	BSKH60373	96T
Human Activin A ELISA Kit	bsk11079	96T
Human ACTN2 ELISA Kit	BSKH60873	96T
Human ACVR1 ELISA Kit	BSKH60874	96T
Human ADAM10 ELISA Kit	BSKH60875	96T
Human ADAM12 ELISA Kit	BSKH60199	96T
Human ADAM17 ELISA Kit	BSKH60247	96T
Human ADAM8 ELISA Kit	BSKH60876	96T
Human ADAMTS5 ELISA Kit	bsk11172	96T
Human Adenosine Deaminase (ADA) ELISA Kit	BSKH61824	96T
Human ADH1 ELISA Kit	BSKH60877	96T
Human Adiponectin ELISA Kit	bsk11015	96T
Human ADIPOR1 ELISA Kit	BSKH60248	96T
Human ADIPOR2 ELISA Kit	BSKH60249	96T
Human ADM ELISA Kit	BSKH60878	96T
Human ADP ELISA Kit	BSKH60181	96T
Human aEP ELISA Kit	BSKH60879	96T
Human AFP ELISA Kit	bsk31008	96T
Human aFP ELISA Kit	BSKH60090	96T
Human AGER ELISA Kit	BSKH60880	96T
Human AGT ELISA Kit	BSKH60250	96T
Human AGTR1 ELISA Kit	BSKH60421	96T
Human AGTR2 ELISA Kit	BSKH60251	96T
Human aHSP ELISA Kit	BSKH60882	96T
Human AIF1 ELISA Kit	BSKH60224	96T
Human ALB ELISA Kit	BSKH60070	96T

产品名称	产品编号	产品规格
Human		
Human ALCAM ELISA Kit	BSKH60252	96T
Human Aldehyde Dehydrogenase, Mitochondrial (ALDM) ELISA Kit	BSKH63087	96T
Human ALDOA ELISA Kit	BSKH60883	96T
Human ALDOS ELISA Kit	BSKH60884	96T
Human Alk-Smase ELISA Kit	BSKH60885	96T
Human ALOX15 ELISA Kit	BSKH60886	96T
Human ALPL ELISA Kit	BSKH61846	96T
Human ALT ELISA Kit	BSKH60887	96T
Human AMH ELISA Kit	BSKH60889	96T
Human Amphiregulin ELISA Kit	bsk11114	96T
Human aMSH ELISA Kit	BSKH60890	96T
Human ANG ELISA Kit	BSKH60145	96T
Human AngI ELISA Kit	BSKH60253	96T
Human AngII ELISA Kit	BSKH60254	96T
Human AngIII ELISA Kit	BSKH60648	96T
Human Angiopoietin-2 ELISA Kit	bsk11115	96T
Human Angiopoietin-like 4 ELISA Kit	bsk11116	96T
Human Angiotensin 1-7 (Ang1-7) ELISA Kit	BSKH62319	96T
Human ANGPT1 ELISA Kit	BSKH60196	96T
Human ANGPT2 ELISA Kit	BSKH60891	96T
Human ANGPT4 ELISA Kit	BSKH60892	96T
Human ANGPTL1 ELISA Kit	BSKH60893	96T
Human Annexin A3 (ANXA3) ELISA Kit	BSKH63103	96T
Human ANP ELISA Kit	BSKH60894	96T
Human Antithrombin (AT) ELISA Kit	BSKH63150	96T
Human ANXA5 ELISA Kit	BSKH60895	96T
Human APOA2 ELISA Kit	BSKH60896	96T
Human APOB ELISA Kit	BSKH60375	96T
Human APOC1 ELISA Kit	BSKH60897	96T
Human APOE ELISA Kit	BSKH60376	96T
Human APOH ELISA Kit	BSKH60898	96T
Human Apolipoprotein A-1/ApoA1 ELISA Kit	bsk11155	96T
Human Apolipoprotein M (APOM) ELISA Kit	BSKH63119	96T
Human AQP1 ELISA Kit	BSKH60899	96T
Human AQP2 ELISA Kit	BSKH60900	96T
Human AQP3 ELISA Kit	BSKH60901	96T
Human AQP4 ELISA Kit	BSKH60902	96T
Human AQP5 ELISA Kit	BSKH60903	96T
Human AQP9 ELISA Kit	BSKH60904	96T
Human Aryl Hydrocarbon Receptor (AhR) ELISA Kit	BSKH61838	96T
Human Asprosin ELISA Kit	BSKH61336	96T
Human ATP1b3 ELISA Kit	BSKH60905	96T
Human ATP4a ELISA Kit	BSKH60649	96T
Human ATXN1 ELISA Kit	BSKH60906	96T
Human b2M ELISA Kit	BSKH60080	96T
Human bACE1 ELISA Kit	BSKH60907	96T

产品名称	产品编号	产品规格
 Human		
Human BAFF/TNFSF13B ELISA Kit	bsk11035	96T
Human BAFFR ELISA Kit	BSKH60908	96T
Human Bcl11A ELISA Kit	BSKH60650	96T
Human Bcl2 ELISA Kit	BSKH60255	96T
Human BDNF ELISA Kit	bsk11023	96T
Human Bid ELISA Kit	BSKH60910	96T
Human Big Dynorphin (Dyn) ELISA Kit	BSKH61976	96T
Human BLK ELISA Kit	BSKH60651	96T
Human BMP1 ELISA Kit	BSKH60911	96T
Human BMP-2 ELISA Kit	bsk11093	96T
Human BMP4 ELISA Kit	BSKH60257	96T
Human BMP-6 ELISA Kit	bsk11127	96T
Human BMP7 ELISA Kit	BSKH60913	96T
Human BMP-7 ELISA Kit	bsk11126	96T
Human BMPR1A ELISA Kit	BSKH60258	96T
Human BMPR1B ELISA Kit	BSKH60914	96T
Human BNP ELISA Kit	BSKH60378	96T
Human bTC ELISA Kit	BSKH60915	96T
Human C Reactive Protein (CRP) ELISA Kit	BSKH60168	96T
Human C1INH ELISA Kit	BSKH60917	96T
Human C1QTNF1 ELISA Kit	BSKH60652	96T
Human C1QTNF3 ELISA Kit	BSKH60653	96T
Human C2 ELISA Kit	BSKH60918	96T
Human C3 ELISA Kit	BSKH60122	96T
Human C4BPb ELISA Kit	BSKH60920	96T
Human C5a ELISA Kit	BSKH60379	96T
Human C7 ELISA Kit	BSKH60921	96T
Human CA1 ELISA Kit	BSKH60923	96T
Human CA-125 ELISA Kit	bsk31012	96T
Human CA15-3 ELISA Kit	BSKH60922	96T
Human CA19-9 ELISA Kit	BSKH60380	96T
Human CA2 ELISA Kit	BSKH60381	96T
Human CA9 ELISA Kit	BSKH60166	96T
Human Calprotectin (CALPRO) ELISA Kit	BSKH63207	96T
Human CAM ELISA Kit	BSKH60925	96T
Human Carbohydrate Antigen 72-4 (CA72-4) ELISA Kit	BSKH60924	96T
Human Carcinoembryonic Antigen (CEA) ELISA Kit	BSKH60091	96T
Human Cartilage Associated Protein (CRTAP) ELISA Kit	BSKH63342	96T
Human CASP14 ELISA Kit	BSKH60926	96T
Human CASP3 ELISA Kit	BSKH60927	96T
Human CASP5 ELISA Kit	BSKH60928	96T
Human CASP7 ELISA Kit	BSKH60382	96T
Human CASP8 ELISA Kit	BSKH60383	96T
Human CASP9 ELISA Kit	BSKH60929	96T
Human CAT ELISA Kit	BSKH60384	96T
Human CAV1 ELISA Kit	BSKH60930	96T
Human CC16 ELISA Kit	BSKH60931	96T

产品名称	产品编号	产品规格
 Human		
Human CCK ELISA Kit	BSKH60933	96T
Human CCKAR ELISA Kit	BSKH60932	96T
Human CCL1 ELISA Kit	BSKH60139	96T
Human CCL1/I-309 ELISA Kit	bsk11150	96T
Human CCL17/TARC ELISA Kit	bsk11053	96T
Human CCL19/MIP-3 beta ELISA Kit	bsk11146	96T
Human CCL21 ELISA Kit	bsk11110	96T
Human CCL22/MDC ELISA Kit	bsk11147	96T
Human CCL24/Eotaxin-2 ELISA Kit	bsk11149	96T
Human CCL25/TECK ELISA Kit	bsk11151	96T
Human CCL26/Eotaxin-3 ELISA Kit	bsk11148	96T
Human CCL3L1 ELISA Kit	BSKH60934	96T
Human CCND1 ELISA Kit	BSKH60935	96T
Human CCND2 ELISA Kit	BSKH60936	96T
Human CCR2 ELISA Kit	BSKH60937	96T
Human CCR3 ELISA Kit	BSKH60939	96T
Human CCR4 ELISA Kit	BSKH60940	96T
Human CCR5 ELISA Kit	BSKH60941	96T
Human CD14 ELISA Kit	BSKH60942	96T
Human CD160 ELISA Kit	BSKH60943	96T
Human CD200 ELISA Kit	BSKH60944	96T
Human CD27/TNFRSF7 ELISA Kit	bsk11167	96T
Human CD276 ELISA Kit	BSKH60945	96T
Human CD30L ELISA Kit	BSKH60946	96T
Human CD5 Antigen Like Protein (CD5L) ELISA Kit	BSKH66009	96T
Human CD8a ELISA Kit	BSKH60385	96T
Human CDH5 ELISA Kit	BSKH60203	96T
Human CDHE ELISA Kit	BSKH60140	96T
Human CDK5 ELISA Kit	BSKH60948	96T
Human CDKN2A ELISA Kit	BSKH60949	96T
Human CEA ELISA Kit	bsk31010	96T
Human CEBPa ELISA Kit	BSKH60950	96T
Human CEBPb ELISA Kit	BSKH60951	96T
Human CENPE ELISA Kit	BSKH60952	96T
Human CENPF ELISA Kit	BSKH60953	96T
Human CETP ELISA Kit	BSKH60954	96T
Human CFD ELISA Kit	BSKH60176	96T
Human CFH ELISA Kit	BSKH60955	96T
Human CFP ELISA Kit	BSKH60956	96T
Human CG ELISA Kit	BSKH60397	96T
Human CGRP ELISA Kit	BSKH60386	96T
Human CHEM ELISA Kit	BSKH60957	96T
Human Cholesterol-25-Hydroxylase (CH25H) ELISA Kit	BSKH63253	96T
Human CHRNa1 ELISA Kit	BSKH60958	96T
Human Chromogranin A (CHGA) ELISA Kit	BSKH61940	96T
Human CKM ELISA Kit	BSKH60959	96T
Human CKMT1A ELISA Kit	BSKH60960	96T

产品名称	产品编号	产品规格
Human		
Human Claudin 5 (CLDN5) ELISA Kit	BSKH63274	96T
Human CLU ELISA Kit	BSKH60204	96T
Human CNP ELISA Kit	BSKH60961	96T
Human CNTF ELISA Kit	BSKH60128	96T
Human CNX ELISA Kit	BSKH60962	96T
Human Coagulation Factor VIII (F8) ELISA Kit	BSKH61994	96T
Human Coenzyme Q10 Homolog B (COQ10B) ELISA Kit	BSKH66027	96T
Human COL1a1 ELISA Kit	BSKH60963	96T
Human COL1a2 ELISA Kit	BSKH60964	96T
Human COL4a1 ELISA Kit	BSKH60965	96T
Human Collagen Type VII (COL7) ELISA Kit	BSKH63299	96T
Human Complement Component 1, R Subcomponent (C1r) ELISA Kit	BSKH61892	96T
Human Complement Component 4a (C4a) ELISA Kit	BSKH60919	96T
Human Complement Component 4b (C4b) ELISA Kit	BSKH61894	96T
Human Complement Component 5 (C5) ELISA Kit	BSKH61896	96T
Human Complement Component C5a ELISA Kit	bsk11176	96T
Human Complement Factor B (CFB) ELISA Kit	BSKH63250	96T
Human Complement Receptor 1, Erythrocyte (CR1) ELISA Kit	BSKH61947	96T
Human Corin ELISA Kit	bsk11121	96T
Human CP ELISA Kit	BSKH60387	96T
Human C-Peptide ELISA Kit	BSKH60217	96T
Human CPP ELISA Kit	BSKH60966	96T
Human CR ELISA Kit	BSKH60969	96T
Human CR2 ELISA Kit	BSKH60967	96T
Human CRH ELISA Kit	BSKH60968	96T
Human CRP ELISA Kit	bsk11041	96T
Human CSF2Rb ELISA Kit	BSKH60422	96T
Human CSH1 ELISA Kit	BSKH60970	96T
Human CST3 ELISA Kit	BSKH60129	96T
Human CSTA ELISA Kit	BSKH60971	96T
Human CT ELISA Kit	BSKH60388	96T
Human CT1 ELISA Kit	BSKH60972	96T
Human CTACK ELISA Kit	BSKH60141	96T
Human CTGF ELISA Kit	BSKH60182	96T
Human CTLA4 ELISA Kit	BSKH60423	96T
Human CTSA ELISA Kit	BSKH60973	96T
Human CTSK ELISA Kit	BSKH60974	96T
Human CTSL ELISA Kit	BSKH60975	96T
Human CTSV ELISA Kit	BSKH60200	96T
Human CTXI ELISA Kit	BSKH60976	96T
Human CTXII ELISA Kit	BSKH60977	96T
Human CTXIII ELISA Kit	BSKH60978	96T
Human Cu/ZnSOD ELISA Kit	bsk11096	96T
Human CX3CL1 ELISA Kit	BSKH60979	96T
Human CX43 ELISA Kit	BSKH60980	96T
Human CXCL1/GRO ELISA Kit	bsk11086	96T

产品名称	产品编号	产品规格
Human		
Human CXCL10/IP10 ELISA Kit	bsk11036	96T
Human CXCL11/I-TAC ELISA Kit	bsk11143	96T
Human CXCL12/SDF-1 ELISA Kit	bsk11087	96T
Human CXCL13/BCA-1/BLC ELISA Kit	bsk11141	96T
Human CXCL14/BRAK ELISA Kit	bsk11145	96T
Human CXCL16 ELISA Kit	bsk11105	96T
Human CXCL2/GRO β ELISA Kit	bsk11164	96T
Human CXCL4/PF4 ELISA Kit	bsk11140	96T
Human CXCL5/ENA-78 ELISA Kit	bsk11052	96T
Human CXCL6/GCP-2 ELISA Kit	bsk11144	96T
Human CXCL7/NAP-2 ELISA Kit	bsk11142	96T
Human CXCR3 ELISA Kit	BSKH60981	96T
Human CXCR4 ELISA Kit	BSKH60982	96T
Human CXCR5 ELISA Kit	BSKH60983	96T
Human CXCR6 ELISA Kit	BSKH60984	96T
Human CYCS ELISA Kit	BSKH61339	96T
Human CYP2E1 ELISA Kit	BSKH60985	96T
Human CYP4 ELISA Kit	BSKH60986	96T
Human CYPB ELISA Kit	BSKH60987	96T
Human D2D ELISA Kit	BSKH60389	96T
Human DAF ELISA Kit	BSKH60988	96T
Human DAO ELISA Kit	BSKH60989	96T
Human DDR1 ELISA Kit	bsk11067	96T
Human DEFB2 ELISA Kit	BSKH60990	96T
Human Deiodinase, Iodothyronine, Type II (DIO2) ELISA Kit	BSKH63419	96T
Human Des ELISA Kit	BSKH60991	96T
Human Diaphanous Homolog 1 (DIAPH1) ELISA Kit	BSKH63416	96T
Human Divalent Metal Transporter 1 (DMT1) ELISA Kit	BSKH63426	96T
Human DKK1 ELISA Kit	BSKH60992	96T
Human DLD ELISA Kit	BSKH60993	96T
Human DNA Methyltransferase 1 (DNMT1) ELISA Kit	BSKH63430	96T
Human DPP4 ELISA Kit	BSKH60994	96T
Human DPPIV/CD26 ELISA Kit	bsk11051	96T
Human DRD2 ELISA Kit	BSKH60995	96T
Human DRD3 ELISA Kit	BSKH60996	96T
Human DSG1 ELISA Kit	BSKH60997	96T
Human DSG3 ELISA Kit	BSKH60998	96T
Human dSIP ELISA Kit	BSKH60999	96T
Human ECE1 ELISA Kit	BSKH61000	96T
Human ECF ELISA Kit	BSKH60160	96T
Human EDN1 ELISA Kit	BSKH60259	96T
Human EDN2 ELISA Kit	BSKH60654	96T
Human EGF ELISA Kit	bsk11025	96T
Human EGFR ELISA Kit	BSKH60260	96T
Human EGFR2 ELISA Kit	BSKH60424	96T
Human EGR1 ELISA Kit	BSKH61001	96T
Human ELA2 ELISA Kit	BSKH61002	96T

产品名称	产品编号	产品规格
 Human		
Human EMILIN1 ELISA Kit	BSKH61003	96T
Human ENG ELISA Kit	BSKH60201	96T
Human eNOS ELISA Kit	bsk11103	96T
Human Epidermal Growth Factor (EGF) ELISA Kit	BSKH60009	96T
Human EPO ELISA Kit	BSKH60390	96T
Human EPOR ELISA Kit	BSKH61004	96T
Human ERb ELISA Kit	BSKH61005	96T
Human ERK2 ELISA Kit	BSKH61006	96T
Human F10 ELISA Kit	BSKH61007	96T
Human F12 ELISA Kit	BSKH61008	96T
Human F2 ELISA Kit	BSKH61009	96T
Human F5 ELISA Kit	BSKH61010	96T
Human F9 ELISA Kit	BSKH61011	96T
Human FABP2 ELISA Kit	BSKH60391	96T
Human FABP6 ELISA Kit	BSKH61012	96T
Human FABP8 ELISA Kit	BSKH61013	96T
Human FAK ELISA Kit	BSKH61014	96T
Human FAS ELISA Kit	BSKH60392	96T
Human Fas Ligand/TNFSF6 ELISA Kit	bsk11132	96T
Human Fas/Apo-1 ELISA Kit	bsk11016	96T
Human FASL ELISA Kit	BSKH60150	96T
Human Fatty Acid Binding Protein 1, Liver (FABP1) ELISA Kit	BSKH61995	96T
Human FBN1 ELISA Kit	BSKH61015	96T
Human FcεRII ELISA Kit	BSKH60218	96T
Human FCN1 ELISA Kit	BSKH61016	96T
Human Ferritin, Light Polypeptide (FTL) ELISA Kit	BSKH63558	96T
Human Fetuin A/AHSG ELISA Kit	bsk11160	96T
Human FGF acidic ELISA Kit	bsk11136	96T
Human FGF basic ELISA Kit	bsk11032	96T
Human FGF1 ELISA Kit	BSKH60047	96T
Human FGF10 ELISA Kit	BSKH60425	96T
Human FGF11 ELISA Kit	BSKH60655	96T
Human FGF12 ELISA Kit	BSKH60656	96T
Human FGF13 ELISA Kit	BSKH60657	96T
Human FGF17 ELISA Kit	BSKH60658	96T
Human FGF18 ELISA Kit	BSKH60659	96T
Human FGF-19 ELISA Kit	bsk11135	96T
Human FGF2 ELISA Kit	BSKH60062	96T
Human FGF-21 ELISA Kit	bsk11134	96T
Human FGF-23 ELISA Kit	bsk11123	96T
Human FGF3 ELISA Kit	BSKH60661	96T
Human FGF4 ELISA Kit	BSKH60262	96T
Human FGF5 ELISA Kit	BSKH60263	96T
Human FGF6 ELISA Kit	BSKH60264	96T
Human FGF7 ELISA Kit	BSKH60057	96T
Human FGF8 ELISA Kit	BSKH60662	96T
Human FGF9 ELISA Kit	BSKH60265	96T

产品名称	产品编号	产品规格
 Human		
Human FGFR1 ELISA Kit	BSKH60426	96T
Human FGFR2 ELISA Kit	BSKH60427	96T
Human FGFR3 ELISA Kit	BSKH60428	96T
Human FGFR4 ELISA Kit	BSKH60429	96T
Human FGL2 ELISA Kit	BSKH61017	96T
Human Fibroblast Activation Protein Alpha (FAPa) ELISA Kit	BSKH63519	96T
Human FIL1d ELISA Kit	BSKH60266	96T
Human Flt3 ELISA Kit	BSKH61018	96T
Human Flt3L ELISA Kit	BSKH60125	96T
Human FN ELISA Kit	BSKH60214	96T
Human Four And A Half LIM Domains Protein 1 (FHL1) ELISA Kit	BSKH63534	96T
Human FOXO1 ELISA Kit	BSKH61019	96T
Human FOXO3 ELISA Kit	BSKH61020	96T
Human FS ELISA Kit	BSKH61021	96T
Human FSH ELISA Kit	bsk31001	96T
Human FUCa1 ELISA Kit	BSKH61022	96T
Human FUT3 ELISA Kit	BSKH61023	96T
Human G6PD ELISA Kit	BSKH61024	96T
Human GaA ELISA Kit	BSKH61025	96T
Human GAL1 ELISA Kit	BSKH61028	96T
Human GAL12 ELISA Kit	BSKH61026	96T
Human GAL14 ELISA Kit	BSKH61027	96T
Human GAL2 ELISA Kit	BSKH61029	96T
Human GAL3 ELISA Kit	BSKH61030	96T
Human GAL4 ELISA Kit	BSKH61031	96T
Human GAL7 ELISA Kit	BSKH61032	96T
Human GAL8 ELISA Kit	BSKH61033	96T
Human GAL9 ELISA Kit	BSKH61034	96T
Human Galectin-3 ELISA Kit	bsk11050	96T
Human GAS6 ELISA Kit	BSKH61035	96T
Human GATA3 ELISA Kit	BSKH61036	96T
Human GCK ELISA Kit	BSKH60393	96T
Human G-CSF ELISA Kit	bsk11017	96T
Human GCSFR ELISA Kit	BSKH60430	96T
Human GDF-15 ELISA Kit	bsk11045	96T
Human GDF5 ELISA Kit	BSKH60394	96T
Human GDF9 ELISA Kit	BSKH61037	96T
Human GDH ELISA Kit	BSKH61038	96T
Human GDNF ELISA Kit	bsk11081	96T
Human GDNF ELISA Kit	BSKH60395	96T
Human GFAP ELISA Kit	BSKH61039	96T
Human GH ELISA Kit	BSKH60152	96T
Human GHRH ELISA Kit	BSKH61040	96T
Human GHRL ELISA Kit	BSKH61041	96T
Human GIP ELISA Kit	BSKH61042	96T
Human GITR ELISA Kit	BSKH60431	96T

产品名称	产品编号	产品规格
 Human		
Human GLa ELISA Kit	BSKH61043	96T
Human GLb ELISA Kit	BSKH61044	96T
Human GLP1 ELISA Kit	BSKH61045	96T
Human Glucagon (GC) ELISA Kit	BSKH62023	96T
Human Glucuronidase Beta (GUSb) ELISA Kit	BSKH62051	96T
Human Glutathione Peroxidase 4 (GPX4) ELISA Kit	BSKH63616	96T
Human GMCSF ELISA Kit	BSKH60011	96T
Human GM-CSF ELISA Kit	bsk11018	96T
Human GnRH ELISA Kit	BSKH61046	96T
Human gp130 ELISA Kit	BSKH60216	96T
Human GPC3 ELISA Kit	BSKH61047	96T
Human GPC4 ELISA Kit	BSKH61048	96T
Human GPI ELISA Kit	BSKH60888	96T
Human GPX1 ELISA Kit	BSKH61049	96T
Human Granzyme A ELISA Kit	bsk11162	96T
Human Granzyme B ELISA Kit	bsk11108	96T
Human Grb10 ELISA Kit	BSKH60663	96T
Human Grb2 ELISA Kit	BSKH60664	96T
Human Growth Regulated Oncogene Gamma (GROg) ELISA Kit	BSKH62046	96T
Human GRP ELISA Kit	BSKH61050	96T
Human GS ELISA Kit	BSKH60396	96T
Human GSK3a ELISA Kit	BSKH61051	96T
Human GSTa1 ELISA Kit	BSKH61052	96T
Human GSTa3 ELISA Kit	BSKH61053	96T
Human GSTa4 ELISA Kit	BSKH61054	96T
Human GSTm1 ELISA Kit	BSKH61055	96T
Human GSTo1 ELISA Kit	BSKH61056	96T
Human GSTt1 ELISA Kit	BSKH61057	96T
Human Guanine Nucleotide-Binding Protein G(i) Subunit Alpha-2 (GNAI2) ELISA Kit	BSKH64528	96T
Human GZMA ELISA Kit	BSKH61058	96T
Human GZMB ELISA Kit	BSKH61059	96T
Human GZMM ELISA Kit	BSKH61060	96T
Human HABP1 ELISA Kit	BSKH61061	96T
Human HAL ELISA Kit	BSKH61062	96T
Human HBEGF ELISA Kit	BSKH60183	96T
Human HB-EGF ELISA Kit	bsk11112	96T
Human HCAM ELISA Kit	BSKH61063	96T
Human HCII ELISA Kit	BSKH61064	96T
Human HDGF ELISA Kit	BSKH60267	96T
Human Hepcidin ELISA Kit	bsk11169	96T
Human HEXa ELISA Kit	BSKH61065	96T
Human Hexokinase 2 (HK2) ELISA Kit	BSKH63658	96T
Human Hexosaminidase B Beta (HEXb) ELISA Kit	BSKH61066	96T
Human HGF ELISA Kit	bsk11034	96T
Human HGFAC ELISA Kit	BSKH60665	96T
Human HGFR ELISA Kit	BSKH60666	96T

产品名称	产品编号	产品规格
 Human		
Human HIF1a ELISA Kit	BSKH61067	96T
Human HK1 ELISA Kit	BSKH61068	96T
Human HMG1 ELISA Kit	BSKH61069	96T
Human HNRPA2B1 ELISA Kit	BSKH61070	96T
Human HO1 ELISA Kit	BSKH61071	96T
Human HO2 ELISA Kit	BSKH61072	96T
Human HPA ELISA Kit	BSKH61073	96T
Human HPRT1 ELISA Kit	BSKH61074	96T
Human Hpt ELISA Kit	BSKH60110	96T
Human HPX ELISA Kit	BSKH60398	96T
Human HSD17B13 ELISA Kit	bsk11182	96T
Human HSF1 ELISA Kit	BSKH61075	96T
Human HSF2 ELISA Kit	BSKH61076	96T
Human HSP40 ELISA Kit	BSKH61077	96T
Human HSP47 ELISA Kit	BSKH61078	96T
Human HSP60 ELISA Kit	bsk11130	96T
Human HSP70/HSPA1A ELISA Kit	bsk11131	96T
Human HSP90b1 ELISA Kit	BSKH61079	96T
Human HSPA8 ELISA Kit	BSKH60399	96T
Human HSPB1 ELISA Kit	BSKH61080	96T
Human HSPb2 ELISA Kit	BSKH61081	96T
Human HSPb6 ELISA Kit	BSKH60400	96T
Human HSPG ELISA Kit	BSKH61083	96T
Human Huntingtin (HTT) ELISA Kit	BSKH63698	96T
Human Hypoxia Inducible Factor 2 Alpha (HIF2a) ELISA Kit	BSKH63657	96T
Human IAPP ELISA Kit	BSKH61084	96T
Human ICAM2 ELISA Kit	BSKH61085	96T
Human ICAM3 ELISA Kit	BSKH61086	96T
Human ICAM5 ELISA Kit	BSKH61087	96T
Human ICOS ELISA Kit	BSKH61088	96T
Human ICOSLG ELISA Kit	BSKH61089	96T
Human IFI30 ELISA Kit	BSKH60667	96T
Human IFNa/bR2 ELISA Kit	BSKH60670	96T
Human IFNa10 ELISA Kit	BSKH60668	96T
Human IFNa2 ELISA Kit	BSKH60268	96T
Human IFNa21 ELISA Kit	BSKH60669	96T
Human IFNa4 ELISA Kit	BSKH60269	96T
Human IFNg ELISA Kit	BSKH60002	96T
Human IFN-α ELISA Kit	bsk11059	96T
Human IFN-β ELISA Kit	bsk11161	96T
Human IFN-γ ELISA Kit	bsk11013	96T
Human IgA ELISA Kit	bsk11098	96T
Human IgE ELISA Kit	bsk11158	96T
Human IGF1 ELISA Kit	BSKH60012	96T
Human IGF1R ELISA Kit	BSKH60432	96T
Human IGF2 ELISA Kit	BSKH60270	96T

产品名称	产品编号	产品规格
 Human		
Human IGF2R ELISA Kit	BSKH60433	96T
Human IGFALS ELISA Kit	BSKH60671	96T
Human IGFBP1 ELISA Kit	BSKH60043	96T
Human IGFBP-1 ELISA Kit	bsk11088	96T
Human IGFBP2 ELISA Kit	BSKH60054	96T
Human IGFBP3 ELISA Kit	BSKH60044	96T
Human IGFBP4 ELISA Kit	BSKH60031	96T
Human IGFBP5 ELISA Kit	BSKH60672	96T
Human IGFBP6 ELISA Kit	BSKH60271	96T
Human IGFBP-7/IGFBP-rp1 ELISA Kit	bsk11157	96T
Human IGF-I ELISA Kit	bsk11058	96T
Human IgG total ELISA Kit	bsk11063	96T
Human IgM ELISA Kit	bsk11064	96T
Human IL-1 R4/ST2 ELISA Kit	bsk11022	96T
Human IL-10 ELISA Kit	bsk11010	96T
Human IL10Ra ELISA Kit	BSKH60435	96T
Human IL-11 ELISA Kit	bsk11071	96T
Human IL11Ra ELISA Kit	BSKH60673	96T
Human IL12A ELISA Kit	BSKH60273	96T
Human IL-12p40 ELISA Kit	bsk11012	96T
Human IL-12p70 ELISA Kit	bsk11011	96T
Human IL12Rb1 ELISA Kit	BSKH60436	96T
Human IL12Rb2 ELISA Kit	BSKH60437	96T
Human IL-13 ELISA Kit	bsk11003	96T
Human IL13Ra2 ELISA Kit	BSKH60438	96T
Human IL-15 ELISA Kit	bsk11029	96T
Human IL16 ELISA Kit	BSKH60275	96T
Human IL-17 ELISA Kit	bsk11028	96T
Human IL17B ELISA Kit	BSKH60439	96T
Human IL-17C ELISA Kit	bsk11054	96T
Human IL17D ELISA Kit	BSKH60675	96T
Human IL-17E ELISA Kit	bsk11055	96T
Human IL17F ELISA Kit	BSKH60440	96T
Human IL17RA ELISA Kit	BSKH60441	96T
Human IL-18 BPa ELISA Kit	bsk11171	96T
Human IL18 ELISA Kit	BSKH60185	96T
Human IL-18 ELISA Kit	bsk11056	96T
Human IL18R1 ELISA Kit	BSKH60442	96T
Human IL18RAP ELISA Kit	BSKH60443	96T
Human IL-19 ELISA Kit	bsk11072	96T
Human IL1e ELISA Kit	BSKH60677	96T
Human IL1h ELISA Kit	BSKH60678	96T
Human IL1R1 ELISA Kit	BSKH60278	96T
Human IL1R2 ELISA Kit	BSKH60279	96T
Human IL1RA ELISA Kit	BSKH60280	96T
Human IL1RAP ELISA Kit	BSKH60679	96T
Human IL1RAPL2 ELISA Kit	BSKH60680	96T

产品名称	产品编号	产品规格
 Human		
Human IL-1 β ELISA Kit	bsk11001	96T
Human IL-1 α ELISA Kit	bsk11039	96T
Human IL2 ELISA Kit	BSKH60075	96T
Human IL-2 ELISA Kit	bsk11002	96T
Human IL20 ELISA Kit	BSKH60063	96T
Human IL-21 ELISA Kit	bsk11083	96T
Human IL21R ELISA Kit	BSKH60681	96T
Human IL-22 ELISA Kit	bsk11030	96T
Human IL22Ra2 ELISA Kit	BSKH60683	96T
Human IL-23 ELISA Kit	bsk11006	96T
Human IL23R ELISA Kit	BSKH60684	96T
Human IL24 ELISA Kit	BSKH60685	96T
Human IL-24/MDA-7 ELISA Kit	bsk11060	96T
Human IL25 ELISA Kit	BSKH60446	96T
Human IL26 ELISA Kit	BSKH60447	96T
Human IL-27 ELISA Kit	bsk11073	96T
Human IL-28/IFN-lambda 2 ELISA Kit	bsk11074	96T
Human IL28A ELISA Kit	BSKH60448	96T
Human IL28B ELISA Kit	BSKH60686	96T
Human IL-29/IFN lambda 1 ELISA Kit	bsk11075	96T
Human IL2Ra ELISA Kit	BSKH60449	96T
Human IL2Rg ELISA Kit	BSKH60450	96T
Human IL-3 ELISA Kit	bsk11020	96T
Human IL-31 ELISA Kit	bsk11076	96T
Human IL-32 α ELISA Kit	bsk11057	96T
Human IL-33 ELISA Kit	bsk11068	96T
Human IL-34 ELISA Kit	bsk11078	96T
Human IL35 ELISA Kit	BSKH60689	96T
Human IL-36 gamma/IL-1F9 ELISA Kit	bsk11179	96T
Human IL-37/IL-1F7 ELISA Kit	bsk11122	96T
Human IL3Ra ELISA Kit	BSKH60454	96T
Human IL-4 ELISA Kit	bsk11004	96T
Human IL-5 ELISA Kit	bsk11005	96T
Human IL5Ra ELISA Kit	BSKH60455	96T
Human IL-6 ELISA Kit	bsk11007	96T
Human IL-6 R alpha ELISA Kit	bsk11137	96T
Human IL-7 ELISA Kit	bsk11111	96T
Human IL7R ELISA Kit	BSKH60457	96T
Human IL-8 ELISA Kit	bsk11008	96T
Human IL8Ra ELISA Kit	BSKH60287	96T
Human IL8Rb ELISA Kit	BSKH60692	96T
Human IL-9 ELISA Kit	bsk11009	96T
Human INHA ELISA Kit	BSKH61090	96T
Human INHbA ELISA Kit	BSKH61091	96T
Human INHbC ELISA Kit	BSKH61092	96T
Human INHbE ELISA Kit	BSKH61093	96T
Human INHBP ELISA Kit	BSKH61094	96T

产品名称	产品编号	产品规格
Human		
Human Inhibitory Subunit Of NF Kappa B Alpha (IkBa) ELISA Kit	BSKH62070	96T
Human INS ELISA Kit	BSKH60066	96T
Human Intercellular Adhesion Molecule 1 (ICAM1) ELISA Kit	BSKH60205	96T
Human Interferon Beta (IFNb) ELISA Kit	BSKH60184	96T
Human Interleukin 1 Alpha (IL1a) ELISA Kit	BSKH60277	96T
Human Interleukin 1 Beta (IL1b) ELISA Kit	BSKH60022	96T
Human iPLA2 ELISA Kit	BSKH61095	96T
Human IRAK1 ELISA Kit	BSKH60458	96T
Human IRAK3 ELISA Kit	BSKH60460	96T
Human IRAK4 ELISA Kit	BSKH60461	96T
Human IRF1 ELISA Kit	BSKH60462	96T
Human IRF3 ELISA Kit	BSKH60463	96T
Human IRF5 ELISA Kit	BSKH60464	96T
Human IRF9 ELISA Kit	BSKH60693	96T
Human ISR ELISA Kit	BSKH61096	96T
Human ITaC ELISA Kit	BSKH60033	96T
Human ITLN1 ELISA Kit	BSKH61097	96T
Human JAM2 ELISA Kit	BSKH61098	96T
Human KGF/FGF-7 ELISA Kit	bsk11152	96T
Human Ki-67/MKI67 ELISA Kit	bsk11133	96T
Human Kim1 ELISA Kit	BSKH61099	96T
Human KLK1 ELISA Kit	BSKH60401	96T
Human KLK10 ELISA Kit	BSKH60189	96T
Human KLK11 ELISA Kit	BSKH61100	96T
Human KLK3 ELISA Kit	BSKH60123	96T
Human KLK4 ELISA Kit	BSKH61101	96T
Human KLK5 ELISA Kit	BSKH61102	96T
Human KLK6 ELISA Kit	BSKH61103	96T
Human KLK7 ELISA Kit	BSKH60402	96T
Human KLK8 ELISA Kit	BSKH61104	96T
Human KLK9 ELISA Kit	BSKH61105	96T
Human Klotho ELISA Kit	bsk11109	96T
Human KRT1 ELISA Kit	BSKH61108	96T
Human KRT14 ELISA Kit	BSKH61106	96T
Human KRT16 ELISA Kit	BSKH61107	96T
Human KRT2 ELISA Kit	BSKH61109	96T
Human KRT4 ELISA Kit	BSKH61110	96T
Human KRT5 ELISA Kit	BSKH61111	96T
Human KRT7 ELISA Kit	BSKH61112	96T
Human KRT9 ELISA Kit	BSKH61113	96T
Human L1CAM ELISA Kit	BSKH61114	96T
Human LAB7-2 ELISA Kit	BSKH60289	96T
Human LAG3 ELISA Kit	BSKH60465	96T
Human LAIR2 ELISA Kit	BSKH61115	96T
Human LAMA1 ELISA Kit	BSKH61125	96T
Human LAMB1 ELISA Kit	BSKH61116	96T

产品名称	产品编号	产品规格
Human		
Human LAP ELISA Kit	BSKH61117	96T
Human LAT ELISA Kit	BSKH61118	96T
Human LCP1 ELISA Kit	BSKH60694	96T
Human LDL ELISA Kit	BSKH60082	96T
Human Lectin Like Oxidized Low Density Lipoprotein Receptor 1 (LOX1) ELISA Kit	BSKH62113	96T
Human Leptin ELISA Kit	bsk11027	96T
Human LFA1a ELISA Kit	BSKH60466	96T
Human LFA3 ELISA Kit	BSKH60290	96T
Human LH ELISA Kit	bsk31002	96T
Human LIF ELISA Kit	bsk11095	96T
Human LIFR ELISA Kit	BSKH61120	96T
Human LIPC ELISA Kit	BSKH61121	96T
Human LIPD ELISA Kit	BSKH61122	96T
Human LIPG ELISA Kit	BSKH61123	96T
Human lipocalin-2/NGAL ELISA Kit	bsk11048	96T
Human LMP2 ELISA Kit	BSKH61124	96T
Human Lpa ELISA Kit	BSKH61126	96T
Human LPHN3 ELISA Kit	BSKH61127	96T
Human LPO ELISA Kit	BSKH61128	96T
Human LpPLA2 ELISA Kit	BSKH61129	96T
Human LTb ELISA Kit	BSKH60467	96T
Human LTBP2 ELISA Kit	BSKH60468	96T
Human LTbR ELISA Kit	BSKH60291	96T
Human LTF ELISA Kit	BSKH60222	96T
Human LY9 ELISA Kit	BSKH60469	96T
Human LY96 ELISA Kit	BSKH60695	96T
Human Lysophosphatidylcholine Acyltransferase 3 (LPCAT3) ELISA Kit	BSKH63803	96T
Human MAG ELISA Kit	BSKH61130	96T
Human MASP2 ELISA Kit	BSKH61131	96T
Human MBP ELISA Kit	BSKH61132	96T
Human MCP-1 ELISA Kit	bsk11040	96T
Human MCP-2/CCL8 ELISA Kit	bsk11085	96T
Human MCP3 ELISA Kit	BSKH60116	96T
Human MCP-3 ELISA Kit	bsk11138	96T
Human MCP4 ELISA Kit	BSKH61134	96T
Human M-CSF ELISA Kit	bsk11080	96T
Human Melatonin Receptor 1A (MTNR1A) ELISA Kit	BSKH63874	96T
Human MEP1a ELISA Kit	BSKH61135	96T
Human MEP1b ELISA Kit	BSKH61136	96T
Human MFGE8 ELISA Kit	BSKH60470	96T
Human MGA ELISA Kit	BSKH61137	96T
Human MGST2 ELISA Kit	BSKH61138	96T
Human MIF ELISA Kit	bsk11047	96T
Human MIG/CXCL9 ELISA Kit	bsk11033	96T
Human MIP1a ELISA Kit	BSKH60126	96T

产品名称	产品编号	产品规格
Human		
Human MIP1b ELISA Kit	BSKH61140	96T
Human MIP3b ELISA Kit	BSKH61141	96T
Human MIP-3α/CCL20 ELISA Kit	bsk11084	96T
Human MIP5 ELISA Kit	BSKH61142	96T
Human MK ELISA Kit	BSKH61143	96T
Human MMP1 ELISA Kit	BSKH61149	96T
Human MMP10 ELISA Kit	BSKH61144	96T
Human MMP11 ELISA Kit	BSKH61145	96T
Human MMP12 ELISA Kit	BSKH61146	96T
Human MMP13 ELISA Kit	BSKH61147	96T
Human MMP-13 ELISA Kit	bsk11125	96T
Human MMP15 ELISA Kit	BSKH61148	96T
Human MMP-2 ELISA Kit	bsk11038	96T
Human MMP23B ELISA Kit	BSKH61150	96T
Human MMP-3 ELISA Kit	bsk11124	96T
Human MMP7 ELISA Kit	BSKH60164	96T
Human MMP8 ELISA Kit	BSKH60112	96T
Human MMP-9 ELISA Kit	bsk11042	96T
Human MOG ELISA Kit	BSKH61151	96T
Human MPIF2 ELISA Kit	BSKH60225	96T
Human MPP2 ELISA Kit	BSKH61152	96T
Human MPP6 ELISA Kit	BSKH61153	96T
Human MPZ ELISA Kit	BSKH61154	96T
Human MR ELISA Kit	BSKH61156	96T
Human MRC2 ELISA Kit	BSKH61155	96T
Human MT3 ELISA Kit	BSKH61157	96T
Human MTL ELISA Kit	BSKH61158	96T
Human MUC1 ELISA Kit	BSKH60403	96T
Human MUC17 ELISA Kit	BSKH61159	96T
Human MUC2 ELISA Kit	BSKH61161	96T
Human MUC20 ELISA Kit	BSKH61160	96T
Human MUC5AC ELISA Kit	BSKH61162	96T
Human MUC5B ELISA Kit	BSKH61163	96T
Human MyD88 ELISA Kit	bsk11174	96T
Human Myeloperoxidase (MPO) ELISA Kit	BSKH60226	96T
Human Myeloperoxidase ELISA Kit	bsk11102	96T
Human MYH1 ELISA Kit	BSKH61164	96T
Human MYH2 ELISA Kit	BSKH61165	96T
Human MYO ELISA Kit	BSKH60067	96T
Human NAA10 ELISA Kit	BSKH61166	96T
Human NAGase ELISA Kit	BSKH61167	96T
Human NAP2 ELISA Kit	BSKH60138	96T
Human NAP3 ELISA Kit	BSKH60133	96T
Human Nerve Growth Factor (NGF) ELISA Kit	BSKH60191	96T
Human NES ELISA Kit	BSKH61168	96T
Human Neurofilament, Light Polypeptide (NEFL) ELISA Kit	BSKH63918	96T

产品名称	产品编号	产品规格
Human		
Human Neutrophil Gelatinase Associated Lipocalin (NGAL) ELISA Kit	BSKH60404	96T
Human NGAL ELISA Kit	bsk51001	96T
Human NGB ELISA Kit	BSKH61169	96T
Human Nidogen (NID) ELISA Kit	BSKH63934	96T
Human NISCH ELISA Kit	BSKH61170	96T
Human NMB ELISA Kit	BSKH61172	96T
Human NOS1 ELISA Kit	BSKH61173	96T
Human NOS2 ELISA Kit	BSKH61174	96T
Human NOS3 ELISA Kit	BSKH61175	96T
Human NOSIP ELISA Kit	BSKH61176	96T
Human NOSTRIN ELISA Kit	BSKH61177	96T
Human NOX1 ELISA Kit	BSKH61178	96T
Human NPHN ELISA Kit	BSKH61179	96T
Human NPS ELISA Kit	BSKH60405	96T
Human NPTX1 ELISA Kit	BSKH61180	96T
Human NPTX2 ELISA Kit	BSKH61181	96T
Human NPY ELISA Kit	BSKH61182	96T
Human NRG1 ELISA Kit	BSKH61183	96T
Human NRP1 ELISA Kit	BSKH61184	96T
Human NSE ELISA Kit	BSKH61185	96T
Human NT3 ELISA Kit	BSKH61186	96T
Human NT4 ELISA Kit	BSKH61187	96T
Human NT-ProANP ELISA Kit	BSKH61188	96T
Human NT-ProBNP ELISA Kit	BSKH60227	96T
Human NTX1 ELISA Kit	BSKH61189	96T
Human OC ELISA Kit	BSKH61190	96T
Human OLFM4 ELISA Kit	BSKH61191	96T
Human ON ELISA Kit	BSKH61192	96T
Human OPG ELISA Kit	bsk11091	96T
Human OPN ELISA Kit	bsk11049	96T
Human OSM ELISA Kit	BSKH61193	96T
Human OSM-R ELISA Kit	bsk11062	96T
Human OX ELISA Kit	BSKH61194	96T
Human OxLDL ELISA Kit	BSKH61195	96T
Human PAI2 ELISA Kit	BSKH61196	96T
Human PAPP A ELISA Kit	BSKH61197	96T
Human PAR2 ELISA Kit	BSKH61198	96T
Human PAR4 ELISA Kit	BSKH61199	96T
Human PARC ELISA Kit	BSKH60144	96T
Human Park7/DJ-1 ELISA Kit	bsk11166	96T
Human PARP ELISA Kit	BSKH61200	96T
Human PBEF/Visfatin ELISA Kit	bsk11061	96T
Human PCK1 ELISA Kit	BSKH61201	96T
Human PCNA ELISA Kit	BSKH61202	96T
Human PCT ELISA Kit	BSKH60407	96T
Human PCX ELISA Kit	BSKH61203	96T

产品名称	产品编号	产品规格
Human		
Human PD-1 ELISA Kit	bsk11139	96T
Human PDCD1LG1 ELISA Kit	BSKH61205	96T
Human PDCD1LG2 ELISA Kit	BSKH61206	96T
Human PDCN ELISA Kit	BSKH61207	96T
Human PDGFA ELISA Kit	BSKH60292	96T
Human PDGF-AA ELISA Kit	bsk11153	96T
Human PDGF-AB ELISA Kit	bsk11089	96T
Human PDGFB ELISA Kit	BSKH60696	96T
Human PDGF-BB ELISA Kit	bsk11090	96T
Human PDGFC ELISA Kit	BSKH60697	96T
Human PDGF-DD ELISA Kit	bsk11154	96T
Human PDGFRa ELISA Kit	BSKH60699	96T
Human PDGFRb ELISA Kit	BSKH60700	96T
Human PDGFRL ELISA Kit	BSKH60701	96T
Human PDHa ELISA Kit	BSKH61208	96T
Human PDP ELISA Kit	BSKH61209	96T
Human PECAM1 ELISA Kit	BSKH61210	96T
Human Perforin 1 (PRF1) ELISA Kit	BSKH62215	96T
Human Periostin/OSF-2 ELISA Kit	bsk11092	96T
Human Peroxiredoxin 2 ELISA Kit	bsk11163	96T
Human PEX2 ELISA Kit	BSKH61211	96T
Human PGA ELISA Kit	BSKH61212	96T
Human PGC ELISA Kit	BSKH61213	96T
Human PGD2S ELISA Kit	BSKH61214	96T
Human PHB ELISA Kit	BSKH61215	96T
Human Phospholipase D (PLD) ELISA Kit	BSKH64073	96T
Human PICP ELISA Kit	BSKH61216	96T
Human PIICP ELISA Kit	BSKH61218	96T
Human PIICP ELISA Kit	BSKH61219	96T
Human PIIINP ELISA Kit	BSKH61220	96T
Human PIINP ELISA Kit	BSKH61221	96T
Human PIK3C2a ELISA Kit	BSKH61222	96T
Human PINP ELISA Kit	BSKH61223	96T
Human PK ELISA Kit	BSKH61228	96T
Human PK2 ELISA Kit	BSKH61224	96T
Human PKBg ELISA Kit	BSKH61225	96T
Human PKCd ELISA Kit	BSKH61226	96T
Human PKCe ELISA Kit	BSKH61227	96T
Human PKM2 ELISA Kit	BSKH61229	96T
Human PKN1 ELISA Kit	BSKH61230	96T
Human PKR ELISA Kit	BSKH61233	96T
Human PKR1 ELISA Kit	BSKH61231	96T
Human PKR2 ELISA Kit	BSKH61232	96T
Human PLAA ELISA Kit	BSKH61234	96T
Human Plasmalemma Vesicle Associated Protein (PLVAP) ELISA Kit	BSKH64086	96T

产品名称	产品编号	产品规格
Human		
Human PLCg1 ELISA Kit	BSKH61235	96T
Human Plg ELISA Kit	BSKH60409	96T
Human PLGF ELISA Kit	bsk11094	96T
Human PLP1 ELISA Kit	BSKH61236	96T
Human Podoplanin (PDPN) ELISA Kit	BSKH64046	96T
Human PON1 ELISA Kit	BSKH61237	96T
Human PPARa ELISA Kit	BSKH61238	96T
Human PPARg ELISA Kit	BSKH61239	96T
Human presepsin ELISA Kit	BSKH61240	96T
Human PRKAa1 ELISA Kit	BSKH61241	96T
Human PRL ELISA Kit	bsk31003	96T
Human PROC ELISA Kit	BSKH61242	96T
Human Pro-Collagen I alpha 1 ELISA Kit	bsk11046	96T
Human PROCR ELISA Kit	BSKH61243	96T
Human Proinsulin ELISA Kit	bsk11077	96T
Human Prolactin ELISA Kit	bsk11165	96T
Human Protein Tyrosine Phosphatase Receptor Type G (PTPRG) ELISA Kit	BSKH64149	96T
Human Protein Tyrosine Phosphatase, Non Receptor Type 1 (PTPN1) ELISA Kit	BSKH64147	96T
Human PROZ ELISA Kit	BSKH61244	96T
Human PRSS1 ELISA Kit	BSKH61245	96T
Human PSA ELISA Kit	bsk31011	96T
Human PSA ELISA Kit	BSKV00012	96T
Human P-selectin ELISA Kit	bsk11043	96T
Human PSMD3 ELISA Kit	BSKH61246	96T
Human PSMD4 ELISA Kit	BSKH61338	96T
Human PTEN ELISA Kit	BSKH60410	96T
Human PTGES ELISA Kit	BSKH61247	96T
Human PTGIS ELISA Kit	BSKH61248	96T
Human PTGS1 ELISA Kit	BSKH61249	96T
Human PTGS2 ELISA Kit	BSKH61250	96T
Human PTH ELISA Kit	BSKH61251	96T
Human PTHR2 ELISA Kit	BSKH61252	96T
Human PTHrP ELISA Kit	BSKH61253	96T
Human PTN ELISA Kit	BSKH60471	96T
Human Puromycin Sensitive Aminopeptidase (PSA) ELISA Kit	BSKH64127	96T
Human PYGL ELISA Kit	BSKH61254	96T
Human PYGM ELISA Kit	BSKH61255	96T
Human RAGE ELISA Kit	BSKH60101	96T
Human RANKL ELISA Kit	BSKH61256	96T
Human Rantes ELISA Kit	bsk11097	96T
Human RAR Related Orphan Receptor Gamma (RORg) ELISA Kit	BSKH64199	96T
Human RARa ELISA Kit	BSKH61257	96T
Human RBP1 ELISA Kit	BSKH61258	96T
Human RBP3 ELISA Kit	BSKH61259	96T
Human RBP4 ELISA Kit	bsk11168	96T

产品名称	产品编号	产品规格
Human		
Human Reelin (RL) ELISA Kit	BSKH64189	96T
Human RELA ELISA Kit	BSKH61260	96T
Human REN ELISA Kit	BSKH60220	96T
Human Resistin ELISA Kit	bsk11069	96T
Human RNASE1 ELISA Kit	BSKH61261	96T
Human RNASEH ELISA Kit	BSKH61262	96T
Human RNASET2 ELISA Kit	BSKH61263	96T
Human S100 Calcium Binding Protein A7 (S100A7) ELISA Kit	BSKH64215	96T
Human S100 Calcium Binding Protein A9 (S100A9) ELISA Kit	BSKH62239	96T
Human S100A1 ELISA Kit	BSKH61265	96T
Human S100A10 ELISA Kit	BSKH60412	96T
Human S100A11 ELISA Kit	BSKH60413	96T
Human S100A12 ELISA Kit	BSKH60414	96T
Human S100A6 ELISA Kit	BSKH60415	96T
Human S100A8 ELISA Kit	BSKH60416	96T
Human S100B ELISA Kit	bsk11066	96T
Human S100B ELISA Kit	BSKV0011	96T
Human S100P ELISA Kit	BSKH60417	96T
Human SAA ELISA Kit	BSKH61266	96T
Human SCCA2 ELISA Kit	BSKH61267	96T
Human sCD40L ELISA Kit	bsk11065	96T
Human SCF ELISA Kit	BSKH60155	96T
Human SCFR ELISA Kit	BSKH60165	96T
Human Secreted Frizzled Related Protein 5 (SFRP5) ELISA Kit	BSKH64250	96T
Human Selectin, Endothelium (SELE) ELISA Kit	BSKH60146	96T
Human Selectin, Platelet (SELP) ELISA Kit	BSKH60127	96T
Human SELL ELISA Kit	BSKH61268	96T
Human Semenogelin I (SEMG1) ELISA Kit	BSKH64236	96T
Human Serine Peptidase Inhibitor Kazal Type 1 (SPINK1) ELISA Kit	BSKH64309	96T
Human Serine/Threonine Kinase 39 (STK39) ELISA Kit	BSKH66205	96T
Human Serpin E1/PAI-1 ELISA Kit	bsk11106	96T
Human sE-selectin ELISA Kit	bsk11044	96T
Human Sestrin 1 (SESN1) ELISA Kit	BSKH64243	96T
Human Sestrin 2 (SESN2) ELISA Kit	BSKH64244	96T
Human Sestrin 3 (SESN3) ELISA Kit	BSKH64245	96T
Human sFRP-5 ELISA Kit	bsk11175	96T
Human SHBG ELISA Kit	BSKH61269	96T
Human sICAM-1 ELISA Kit	bsk11026	96T
Human SIGLEC1 ELISA Kit	BSKH61270	96T
Human Signal Transducer And Activator Of Transcription 1 (STAT1) ELISA Kit	BSKH62263	96T
Human Signal Transducer And Activator Of Transcription 5A (STAT5A) ELISA Kit	BSKH62268	96T
Human Signal Transducer And Activator Of Transcription 5B (STAT5B) ELISA Kit	BSKH62269	96T

产品名称	产品编号	产品规格
Human		
Human Signal Transducer And Activator Of Transcription 6 (STAT6) ELISA Kit	BSKH62270	96T
Human sIL-2R ELISA Kit	bsk11019	96T
Human SIRPa ELISA Kit	BSKH61271	96T
Human SIRPb1 ELISA Kit	BSKH61272	96T
Human SIRPg ELISA Kit	BSKH61273	96T
Human SIRT2 ELISA Kit	BSKH61274	96T
Human SLAMF7 ELISA Kit	BSKH60702	96T
Human Slit2 ELISA Kit	BSKH61276	96T
Human Smad1 ELISA Kit	BSKH61277	96T
Human Smad7 ELISA Kit	BSKH61278	96T
Human SNCg ELISA Kit	BSKH61279	96T
Human SOD3 ELISA Kit	BSKH61280	96T
Human SOD4 ELISA Kit	BSKH61281	96T
Human Sodium Iodide Symporter (NIS) ELISA Kit	BSKH63936	96T
Human Solute Carrier Family 26, Member 4 (SLC26A4) ELISA Kit	BSKH64276	96T
Human Solute Carrier Family 3, Member 2 (SLC3A2) ELISA Kit	BSKH64285	96T
Human SOX1 ELISA Kit	BSKH61283	96T
Human SOX18 ELISA Kit	BSKH61282	96T
Human SP ELISA Kit	BSKH61340	96T
Human SPA ELISA Kit	BSKH61284	96T
Human SPARCL1 ELISA Kit	BSKH60198	96T
Human SPD ELISA Kit	BSKH60207	96T
Human SPTAN1 ELISA Kit	BSKH61285	96T
Human Squalene Epoxidase (SQLE) ELISA Kit	BSKH64318	96T
Human Squamous Cell Carcinoma Antigen 1 (SCCA1) ELISA Kit	BSKH62244	96T
Human SST ELISA Kit	BSKH61286	96T
Human Stanniocalcin 1 (STC1) ELISA Kit	BSKH64330	96T
Human sTNF RI/TNFRSF1A ELISA Kit	bsk11119	96T
Human sTNF RII/TNFRSF1B ELISA Kit	bsk11118	96T
Human Superoxide Dismutase 1, Soluble (SOD1) ELISA Kit	BSKH62259	96T
Human Suppressors Of Cytokine Signaling 1 (SOCS1) ELISA Kit	BSKH64302	96T
Human Suppressors Of Cytokine Signaling 3 (SOCS3) ELISA Kit	BSKH62258	96T
Human Surv ELISA Kit	BSKH60192	96T
Human Syndecan-1 ELISA Kit	bsk11178	96T
Human SYP ELISA Kit	BSKH61287	96T
Human T3 ELISA Kit	bsk31006	96T
Human T4 ELISA Kit	bsk31007	96T
Human TAC1 ELISA Kit	BSKH61171	96T
Human TACR1 ELISA Kit	BSKH61288	96T
Human TAFI ELISA Kit	BSKH61289	96T
Human TAP ELISA Kit	BSKH61290	96T
Human Tar DNA Binding Protein 43kDa (TDP43) ELISA Kit	BSKH62275	96T

产品名称	产品编号	产品规格
Human		
Human Target Of The Antiproliferative Antibody 1 (TAPA1) ELISA Kit	BSKH62272	96T
Human TBG ELISA Kit	BSKH61292	96T
Human TDGF1 ELISA Kit	BSKH60703	96T
Human TEP1 ELISA Kit	BSKH61293	96T
Human TESTO ELISA Kit	bsk31004	96T
Human Tetraspanin 30 (TSPAN30) ELISA Kit	BSKH62296	96T
Human TF ELISA Kit	BSKH61295	96T
Human TFF2 ELISA Kit	BSKH61294	96T
Human TFPI ELISA Kit	BSKH60147	96T
Human TFR2 ELISA Kit	BSKH61296	96T
Human TG ELISA Kit	BSKH61297	96T
Human TGFb2 ELISA Kit	BSKH60295	96T
Human TGFb1 ELISA Kit	BSKH60296	96T
Human TGFbR1 ELISA Kit	BSKH60297	96T
Human TGFbR2 ELISA Kit	BSKH60704	96T
Human TGF-α ELISA Kit	bsk11104	96T
Human TGF-β1 ELISA Kit	bsk11021	96T
Human TGF-β3 ELISA Kit	bsk11120	96T
Human THBS1 ELISA Kit	BSKH61298	96T
Human THBS2 ELISA Kit	BSKH60180	96T
Human Thrombomodulin/BDCA-3 ELISA Kit	bsk11113	96T
Human Thrombopoietin ELISA Kit	bsk11107	96T
Human Thrombospondin-1 ELISA Kit	bsk11159	96T
Human Tie1 ELISA Kit	BSKH60136	96T
Human Tie2 ELISA Kit	BSKH60175	96T
Human TIEG1 ELISA Kit	BSKH61299	96T
Human TIM-1/KIM-1 ELISA Kit	bsk11070	96T
Human TIMP-1 ELISA Kit	bsk11100	96T
Human TIMP-2 ELISA Kit	bsk11099	96T
Human TIMP3 ELISA Kit	BSKH61301	96T
Human TIMP4 ELISA Kit	BSKH60148	96T
Human TL1A/TNFSF15 ELISA Kit	bsk11170	96T
Human TLN ELISA Kit	BSKH61302	96T
Human TLR2 ELISA Kit	bsk11128	96T
Human TLR4 ELISA Kit	BSKH61304	96T
Human TLR6 ELISA Kit	BSKH61305	96T
Human TLR9 ELISA Kit	BSKH61306	96T
Human TM ELISA Kit	BSKH61308	96T
Human TNFaIP2 ELISA Kit	BSKH60705	96T
Human TNFSF12 ELISA Kit	BSKH60706	96T
Human TNFSF13 ELISA Kit	BSKH60474	96T
Human TNFSF14 ELISA Kit	BSKH60299	96T
Human TNFSF4 ELISA Kit	BSKH60300	96T
Human TNFSF7 ELISA Kit	BSKH60475	96T
Human TNFSF9 ELISA Kit	BSKH60707	96T
Human TNF-α ELISA Kit	bsk11014	96T

产品名称	产品编号	产品规格
Human		
Human TNF-β ELISA Kit	bsk11156	96T
Human TNNI3 ELISA Kit	BSKH60418	96T
Human TOP2 ELISA Kit	BSKH61309	96T
Human Total MMP-1 ELISA Kit	bsk11037	96T
Human TP ELISA Kit	BSKH61311	96T
Human TP53 ELISA Kit	BSKH61310	96T
Human tPA ELISA Kit	BSKH60419	96T
Human TPO ELISA Kit	BSKH61314	96T
Human TPPP ELISA Kit	BSKH61313	96T
Human TRAIL ELISA Kit	BSKH60040	96T
Human TRAIL/TNFSF10 ELISA Kit	bsk11173	96T
Human TRANCE/RANK L/TNFSF11 ELISA Kit	bsk11177	96T
Human Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit	BSKH60051	96T
Human TREM1 ELISA Kit	BSKH61315	96T
Human TRF ELISA Kit	BSKH60095	96T
Human TRH ELISA Kit	BSKH61316	96T
Human Tripeptidyl Peptidase I (TPP1) ELISA Kit	BSKH64422	96T
Human Trx ELISA Kit	BSKH61317	96T
Human TrxR1 ELISA Kit	BSKH61318	96T
Human TSC22 ELISA Kit	BSKH60476	96T
Human TSH ELISA Kit	bsk31005	96T
Human TSLP ELISA Kit	bsk11031	96T
Human TTR ELISA Kit	BSKH61319	96T
Human Tubulin Beta (TUBb) ELISA Kit	BSKH62298	96T
Human TYR ELISA Kit	BSKH61320	96T
Human Tyrosine Kinase 2 (Tyk2) ELISA Kit	BSKH62300	96T
Human Ub ELISA Kit	BSKH61321	96T
Human Ubiquitin Carboxyl Terminal Hydrolase L1 (UCHL1) ELISA Kit	BSKH64457	96T
Human UCN1 ELISA Kit	BSKH61322	96T
Human uPA ELISA Kit	BSKH60420	96T
Human uPAR ELISA Kit	BSKH61323	96T
Human UROC1 ELISA Kit	BSKH61324	96T
Human Uromodulin ELISA Kit	bsk11180	96T
Human UST2 ELISA Kit	BSKH61325	96T
Human UTS2R ELISA Kit	BSKH61326	96T
Human Vaccinia Related Kinase 1 (VRK1) ELISA Kit	BSKH64488	96T
Human Vaspin ELISA Kit	BSKH61327	96T
Human VCAM-1 ELISA Kit	bsk11101	96T
Human VDR ELISA Kit	BSKH61328	96T
Human VEGF ELISA Kit	bsk11024	96T
Human VEGF145 ELISA Kit	BSKH60078	96T
Human VEGF165 ELISA Kit	BSKH60477	96T
Human VEGFB ELISA Kit	BSKH60301	96T
Human VEGFC ELISA Kit	BSKH60064	96T
Human VEGFD ELISA Kit	BSKH60052	96T
Human VEGFR1 ELISA Kit	BSKH60065	96T

产品名称	产品编号	产品规格
 Human		
Human VEGFR2 ELISA Kit	BSKH60478	96T
Human VEGFR3 ELISA Kit	BSKH60479	96T
Human VF ELISA Kit	BSKH61329	96T
Human VGF ELISA Kit	BSKH60480	96T
Human VIP ELISA Kit	BSKH61330	96T
Human Vitamin D BP ELISA Kit	bsk11129	96T
Human vWA3A ELISA Kit	BSKH61331	96T
Human vWF ELISA Kit	BSKH61332	96T
Human vWF-A2 ELISA Kit	bsk11117	96T
Human vWFCP ELISA Kit	BSKH60221	96T
Human WFDC2 ELISA Kit	BSKH61333	96T
Human Wilms Tumor Protein (WT1) ELISA Kit	BSKH64495	96T
Human XPC ELISA Kit	BSKH61334	96T
Human XRCC6 ELISA Kit	BSKH61335	96T
Human β -Thromboglobulin (β TG) ELISA Kit	BSKH60916	96T
Human β -HCG ELISA Kit	bsk31009	96T
Human β -NGF ELISA Kit	bsk11082	96T
 Mouse		
Mouse a1AGP ELISA Kit	BSKM61350	96T
Mouse a1M ELISA Kit	BSKM61351	96T
Mouse a2M ELISA Kit	BSKM60232	96T
Mouse a2PI ELISA Kit	BSKM61352	96T
Mouse ABCG2 ELISA Kit	BSKM61355	96T
Mouse ACE ELISA Kit	BSKM60484	96T
Mouse Acetylcholinesterase (AChE) ELISA Kit	BSKM62369	96T
Mouse ACP5 ELISA Kit	BSKM61356	96T
Mouse ACTH ELISA Kit	BSKM60485	96T
Mouse ADAM8 ELISA Kit	BSKM61357	96T
Mouse ADH1 ELISA Kit	BSKM61358	96T
Mouse ADIPOR1 ELISA Kit	BSKM60311	96T
Mouse ADM ELISA Kit	BSKM61359	96T
Mouse ADP ELISA Kit	BSKM60008	96T
Mouse AGER ELISA Kit	BSKM61362	96T
Mouse AGT ELISA Kit	BSKM60312	96T
Mouse aHSG ELISA Kit	BSKM61363	96T
Mouse ALB ELISA Kit	BSKM60042	96T
Mouse ALCAM ELISA Kit	BSKM60313	96T
Mouse ALDOA ELISA Kit	BSKM61364	96T
Mouse ALDOS ELISA Kit	BSKM61365	96T
Mouse Alk-Smase ELISA Kit	BSKM61366	96T
Mouse Alpha-Fetoprotein (aFP) ELISA Kit	BSKM61361	96T
Mouse ALT ELISA Kit	BSKM61367	96T
Mouse AMH ELISA Kit	BSKM61369	96T
Mouse AngI ELISA Kit	BSKM60315	96T
Mouse AngII ELISA Kit	BSKM60314	96T

产品名称	产品编号	产品规格
 Mouse		
Mouse Angiotensin I Converting Enzyme 2 (ACE2) ELISA Kit	BSKM62366	96T
Mouse ANGPT1 ELISA Kit	BSKM61370	96T
Mouse ANGPT2 ELISA Kit	BSKM61371	96T
Mouse ANGPTL3 ELISA Kit	BSKM60100	96T
Mouse ANP ELISA Kit	BSKM61372	96T
Mouse ANXA5 ELISA Kit	BSKM61373	96T
Mouse APOA1 ELISA Kit	BSKM60486	96T
Mouse APOA2 ELISA Kit	BSKM61374	96T
Mouse APOC1 ELISA Kit	BSKM61375	96T
Mouse APOE ELISA Kit	BSKM61376	96T
Mouse APOH ELISA Kit	BSKM61377	96T
Mouse AQP1 ELISA Kit	BSKM61378	96T
Mouse AQP2 ELISA Kit	BSKM61379	96T
Mouse AQP3 ELISA Kit	BSKM61380	96T
Mouse AQP4 ELISA Kit	BSKM61381	96T
Mouse AQP5 ELISA Kit	BSKM61382	96T
Mouse AREG ELISA Kit	BSKM61383	96T
Mouse Aspartate Aminotransferase (AST) ELISA Kit	BSKM62425	96T
Mouse ATP1b3 ELISA Kit	BSKM61384	96T
Mouse ATP4a ELISA Kit	BSKM60769	96T
Mouse b2M ELISA Kit	BSKM61385	96T
Mouse bACE1 ELISA Kit	BSKM61386	96T
Mouse BAFF/BLyS/TNFSF13B ELISA Kit	bsk12049	96T
Mouse BAFFR ELISA Kit	BSKM61387	96T
Mouse Bcl2 ELISA Kit	BSKM60316	96T
Mouse Bcl2L ELISA Kit	BSKM60770	96T
Mouse Bcl9 ELISA Kit	BSKM60771	96T
Mouse Beta-Thromboglobulin (β TG) ELISA Kit	BSKM61393	96T
Mouse Bid ELISA Kit	BSKM60487	96T
Mouse Big Dynorphin (Dyn) ELISA Kit	BSKM61977	96T
Mouse BLC1 ELISA Kit	BSKM60018	96T
Mouse BMP1 ELISA Kit	BSKM61389	96T
Mouse BMP2 ELISA Kit	BSKM60318	96T
Mouse BMP4 ELISA Kit	BSKM60319	96T
Mouse BMP6 ELISA Kit	BSKM61390	96T
Mouse BMP7 ELISA Kit	BSKM61391	96T
Mouse BMPR1A ELISA Kit	BSKM60320	96T
Mouse BNP ELISA Kit	BSKM61392	96T
Mouse Brain Derived Neurotrophic Factor (BDNF) ELISA Kit	BSKM60317	96T
Mouse bTC ELISA Kit	BSKM60307	96T
Mouse C1INH ELISA Kit	BSKM61394	96T
Mouse C1qA ELISA Kit	BSKM66484	96T
Mouse C3 ELISA Kit	BSKM61395	96T
Mouse C4a ELISA Kit	BSKM61396	96T
Mouse C5a ELISA Kit	BSKM60488	96T
Mouse CA1 ELISA Kit	BSKM61398	96T
Mouse Calcitonin Gene Related Peptide (CGRP) ELISA Kit	BSKM60491	96T

产品名称	产品编号	产品规格
 Mouse		
Mouse Calreticulin (CRT) ELISA Kit	BSKM62523	96T
Mouse Carbohydrate Antigen 125 (CA125) ELISA Kit	BSKM61397	96T
Mouse CART ELISA Kit	BSKM61399	96T
Mouse CASP14 ELISA Kit	BSKM61400	96T
Mouse CASP3 ELISA Kit	BSKM61401	96T
Mouse CASP7 ELISA Kit	BSKM60489	96T
Mouse CASP8 ELISA Kit	BSKM60490	96T
Mouse CASP9 ELISA Kit	BSKM61402	96T
Mouse CAV1 ELISA Kit	BSKM61403	96T
Mouse CC16 ELISA Kit	BSKM61404	96T
Mouse CCKAR ELISA Kit	BSKM61405	96T
Mouse CCL11/Eotaxin-1 ELISA Kit	bsk12063	96T
Mouse CCL24/Eotaxin/MPIF-2 ELISA Kit	bsk12058	96T
Mouse CCL5/Rantes ELISA Kit	bsk12043	96T
Mouse CCND1 ELISA Kit	BSKM61407	96T
Mouse CCR2 ELISA Kit	BSKM60938	96T
Mouse CD14 ELISA Kit	BSKM61408	96T
Mouse CD30L ELISA Kit	BSKM61409	96T
Mouse CD40L ELISA Kit	BSKM61410	96T
Mouse CD8a ELISA Kit	BSKM66479	96T
Mouse CDHE ELISA Kit	BSKM61411	96T
Mouse CDKN2A ELISA Kit	BSKM61412	96T
Mouse CEA ELISA Kit	BSKM60233	96T
Mouse CEBPb ELISA Kit	BSKM61413	96T
Mouse CFH ELISA Kit	BSKM61414	96T
Mouse Checkpoint Homolog (CHEK1) ELISA Kit	BSKM66285	96T
Mouse Chemerin ELISA Kit	bsk12074	96T
Mouse Cholecystokinin (CCK) ELISA Kit	BSKM61406	96T
Mouse CKM ELISA Kit	BSKM61415	96T
Mouse CLU ELISA Kit	BSKM60149	96T
Mouse CNP ELISA Kit	BSKM61416	96T
Mouse CNTF ELISA Kit	BSKM61417	96T
Mouse COL1a1 ELISA Kit	BSKM61418	96T
Mouse COL1a2 ELISA Kit	BSKM61419	96T
Mouse COL4a1 ELISA Kit	BSKM61420	96T
Mouse COMP ELISA Kit	BSKM60492	96T
Mouse Complement Component 4b (C4b) ELISA Kit	BSKM62452	96T
Mouse Complement Component C5a ELISA Kit	bsk12057	96T
Mouse CP ELISA Kit	BSKM61421	96T
Mouse C-Peptide ELISA Kit	BSKM60493	96T
Mouse CRP ELISA Kit	bsk12069	96T
Mouse CST3 ELISA Kit	BSKM61425	96T
Mouse CSTA ELISA Kit	BSKM61426	96T
Mouse CT ELISA Kit	BSKM61428	96T
Mouse CT1 ELISA Kit	BSKM61427	96T
Mouse CTGF ELISA Kit	BSKM60321	96T
Mouse CTSA ELISA Kit	BSKM61429	96T

产品名称	产品编号	产品规格
 Mouse		
Mouse CTSK ELISA Kit	BSKM61430	96T
Mouse CTSL ELISA Kit	BSKM61431	96T
Mouse CX3CL1/Fractalkine ELISA Kit	bsk12080	96T
Mouse CX43 ELISA Kit	BSKM61434	96T
Mouse CXCL1 ELISA Kit	BSKM1005	96T
Mouse CXCL10/IP10 ELISA Kit	bsk12045	96T
Mouse CXCL11/I-TAC ELISA Kit	bsk12059	96T
Mouse CXCL12/SDF-1 ELISA Kit	bsk12056	96T
Mouse CXCL15/Lungkine ELISA Kit	bsk12061	96T
Mouse CXCL7/Thymus Chemokine-1 ELISA Kit	bsk12060	96T
Mouse CXCL9/MIG ELISA Kit	bsk12044	96T
Mouse CXCR3 ELISA Kit	BSKM61435	96T
Mouse CXCR4 ELISA Kit	BSKM61436	96T
Mouse Cyclin Dependent Kinase Inhibitor 1A (CDKN1A) ELISA Kit	BSKM64659	96T
Mouse CYCS ELISA Kit	BSKM61437	96T
Mouse CYP2E1 ELISA Kit	BSKM61438	96T
Mouse CYP4 ELISA Kit	BSKM61439	96T
Mouse CYPB ELISA Kit	BSKM61440	96T
Mouse DDR1 ELISA Kit	BSKM61441	96T
Mouse DEFb2 ELISA Kit	BSKM61442	96T
Mouse Des ELISA Kit	BSKM61443	96T
Mouse Diazepam Binding Inhibitor (DBI) ELISA Kit	BSKM64720	96T
Mouse DKK1 ELISA Kit	BSKM60170	96T
Mouse DLAT ELISA Kit	BSKM61444	96T
Mouse DLD ELISA Kit	BSKM61445	96T
Mouse DNA Methyltransferase 1 (DNMT1) ELISA Kit	BSKM64734	96T
Mouse DPP4 ELISA Kit	BSKM61446	96T
Mouse EDN2 ELISA Kit	BSKM66471	96T
Mouse EGF ELISA Kit	bsk12036	96T
Mouse EGFR ELISA Kit	BSKM60323	96T
Mouse EGFR2 ELISA Kit	BSKM60531	96T
Mouse EGR1 ELISA Kit	BSKM61448	96T
Mouse EG-VEGF ELISA Kit	BSKM60324	96T
Mouse ENG ELISA Kit	BSKM60130	96T
Mouse Enolase, Neuron Specific (NSE) ELISA Kit	BSKM61580	96T
Mouse Epiregulin (EREG) ELISA Kit	BSKM62568	96T
Mouse ERb ELISA Kit	BSKM61450	96T
Mouse Erythropoietin ELISA Kit	bsk12081	96T
Mouse E-Selectin/CD62E ELISA Kit	bsk12067	96T
Mouse eutrophil Elastase/ELA2 ELISA Kit	bsk12076	96T
Mouse F10 ELISA Kit	BSKM61451	96T
Mouse F12 ELISA Kit	BSKM61452	96T
Mouse F2 ELISA Kit	BSKM61453	96T
Mouse F5 ELISA Kit	BSKM61454	96T
Mouse F9 ELISA Kit	BSKM61455	96T
Mouse FABP2 ELISA Kit	BSKM60494	96T

产品名称	产品编号	产品规格
 Mouse		
Mouse FABP6 ELISA Kit	BSKM61456	96T
Mouse FAS ELISA Kit	BSKM60495	96T
Mouse FASL ELISA Kit	BSKM60151	96T
Mouse FCN1 ELISA Kit	BSKM61457	96T
Mouse FFAR3 ELISA Kit	BSKM61458	96T
Mouse FGF basic ELISA Kit	bsk12037	96T
Mouse FGF1 ELISA Kit	BSKM60325	96T
Mouse FGF10 ELISA Kit	BSKM60534	96T
Mouse FGF11 ELISA Kit	BSKM60777	96T
Mouse FGF12 ELISA Kit	BSKM60778	96T
Mouse FGF13 ELISA Kit	BSKM60779	96T
Mouse FGF18 ELISA Kit	BSKM60780	96T
Mouse FGF2 ELISA Kit	BSKM60327	96T
Mouse FGF-21 ELISA Kit	bsk12010	96T
Mouse FGF23 ELISA Kit	BSKM60326	96T
Mouse FGF3 ELISA Kit	BSKM60782	96T
Mouse FGF4 ELISA Kit	BSKM60328	96T
Mouse FGF6 ELISA Kit	BSKM60329	96T
Mouse FGF7 ELISA Kit	BSKM60330	96T
Mouse FGF8 ELISA Kit	BSKM60783	96T
Mouse FGF9 ELISA Kit	BSKM60331	96T
Mouse FGFR1 ELISA Kit	BSKM60546	96T
Mouse FGFR3 ELISA Kit	BSKM60547	96T
Mouse FGL2 ELISA Kit	BSKM61459	96T
Mouse Fibrinogen Alpha (FGa) ELISA Kit	BSKM62586	96T
Mouse Fibrinogen Beta (FGb) ELISA Kit	BSKM62587	96T
Mouse Flt3L ELISA Kit	BSKM60171	96T
Mouse FN ELISA Kit	BSKM61460	96T
Mouse FOXO3 ELISA Kit	BSKM61461	96T
Mouse FSH ELISA Kit	BSKM61462	96T
Mouse G6PD ELISA Kit	BSKM61464	96T
Mouse GAL1 ELISA Kit	BSKM61466	96T
Mouse GAL12 ELISA Kit	BSKM61465	96T
Mouse GAL2 ELISA Kit	BSKM61467	96T
Mouse GAL3 ELISA Kit	BSKM61468	96T
Mouse GAL4 ELISA Kit	BSKM61469	96T
Mouse GAL6 ELISA Kit	BSKM61470	96T
Mouse GAL7 ELISA Kit	BSKM61471	96T
Mouse GAL8 ELISA Kit	BSKM61472	96T
Mouse GAL9 ELISA Kit	BSKM61473	96T
Mouse GAS6 ELISA Kit	BSKM61474	96T
Mouse Gastrin (GT) ELISA Kit	BSKM62631	96T
Mouse GCK ELISA Kit	BSKM61475	96T
Mouse G-CSF ELISA Kit	bsk12031	96T
Mouse GDF-15 ELISA Kit	bsk12038	96T
Mouse GDH ELISA Kit	BSKM61476	96T

产品名称	产品编号	产品规格
 Mouse		
Mouse GDNF ELISA Kit	BSKM61477	96T
Mouse GFAP ELISA Kit	BSKM61478	96T
Mouse GFER ELISA Kit	BSKM60784	96T
Mouse GH ELISA Kit	BSKM61479	96T
Mouse GHRL ELISA Kit	BSKM61481	96T
Mouse GIP ELISA Kit	BSKM61482	96T
Mouse Glucagon Like Peptide 1 (GLP1) ELISA Kit	BSKM61483	96T
Mouse Glucuronidase Beta (GUSb) ELISA Kit	BSKM62632	96T
Mouse Glycogen Synthase Kinase 3 Beta (GSK3b) ELISA Kit	BSKM64819	96T
Mouse Glycoprotein Ib Beta Polypeptide, Platelet (GP1Bb) ELISA Kit	BSKM62620	96T
Mouse GM-CSF ELISA Kit	bsk12028	96T
Mouse GnRH ELISA Kit	BSKM61484	96T
Mouse gp130 ELISA Kit	BSKM61485	96T
Mouse GPC4 ELISA Kit	BSKM61486	96T
Mouse GPI ELISA Kit	BSKM61368	96T
Mouse GPX1 ELISA Kit	BSKM61487	96T
Mouse Granzyme B ELISA Kit	bsk12050	96T
Mouse Growth Associated Protein 43 (GAP43) ELISA Kit	BSKM62600	96T
Mouse Growth Regulated Oncogene Gamma (GROg) ELISA Kit	BSKM62626	96T
Mouse GS ELISA Kit	BSKM61489	96T
Mouse GSTa1 ELISA Kit	BSKM61490	96T
Mouse GSTa3 ELISA Kit	BSKM61491	96T
Mouse GSTa4 ELISA Kit	BSKM61492	96T
Mouse GSTo1 ELISA Kit	BSKM61493	96T
Mouse HBEGF ELISA Kit	BSKM60210	96T
Mouse HDGF ELISA Kit	BSKM60333	96T
Mouse Heat Shock 70kDa Protein 1A (HSPA1A) ELISA Kit	BSKM62653	96T
Mouse Hepcidin (Hepc) ELISA Kit	BSKM62640	96T
Mouse HGF ELISA Kit	bsk12046	96T
Mouse High Mobility Group Protein 1 (HMG1) ELISA Kit	BSKM61502	96T
Mouse HSPA8 ELISA Kit	BSKM60496	96T
Mouse HSPb6 ELISA Kit	BSKM60497	96T
Mouse ICAM1 ELISA Kit	BSKM60096	96T
Mouse IFNa ELISA Kit	BSKM60083	96T
Mouse IFN-beta ELISA Kit	bsk12073	96T
Mouse IFN-γ ELISA Kit	bsk12001	96T
Mouse IgA ELISA Kit	bsk12051	96T
Mouse IgE ELISA Kit	bsk12030	96T
Mouse IGF1 ELISA Kit	BSKM60013	96T
Mouse IGF1R ELISA Kit	BSKM66488	96T
Mouse IGF2 ELISA Kit	BSKM60336	96T
Mouse IGFBP1 ELISA Kit	BSKM60337	96T
Mouse IGFBP2 ELISA Kit	BSKM60338	96T
Mouse IGFBP3 ELISA Kit	BSKM60037	96T
Mouse IGFBP4 ELISA Kit	BSKM60339	96T
Mouse IGFBP6 ELISA Kit	BSKM60340	96T

产品名称	产品编号	产品规格
 Mouse		
Mouse IGFBP7 ELISA Kit	BSKM60566	96T
Mouse IGF-I ELISA Kit	bsk12047	96T
Mouse IgG total ELISA Kit	bsk12052	96T
Mouse IgM ELISA Kit	bsk12053	96T
Mouse IL-10 ELISA Kit	bsk12007	96T
Mouse IL10Ra ELISA Kit	BSKM60569	96T
Mouse IL11 ELISA Kit	BSKM60341	96T
Mouse IL11Ra ELISA Kit	BSKM60795	96T
Mouse IL-12/23 (P40) ELISA Kit	bsk12020	96T
Mouse IL12A ELISA Kit	BSKM60342	96T
Mouse IL-12p70 ELISA Kit	bsk12017	96T
Mouse IL-13 ELISA Kit	bsk12006	96T
Mouse IL13Ra2 ELISA Kit	BSKM60577	96T
Mouse IL-15 ELISA Kit	bsk12023	96T
Mouse IL16 ELISA Kit	BSKM60347	96T
Mouse IL17 ELISA Kit	BSKM60348	96T
Mouse IL-17A ELISA Kit	bsk12021	96T
Mouse IL17B ELISA Kit	BSKM60581	96T
Mouse IL17C ELISA Kit	BSKM60798	96T
Mouse IL-17F ELISA Kit	bsk12065	96T
Mouse IL18 ELISA Kit	BSKM60349	96T
Mouse IL18BP ELISA Kit	BSKM60799	96T
Mouse IL19 ELISA Kit	BSKM60586	96T
Mouse IL1R1 ELISA Kit	BSKM60351	96T
Mouse IL1R2 ELISA Kit	BSKM60352	96T
Mouse IL-1ra/IL-1F3 ELISA Kit	bsk12066	96T
Mouse IL1RAP ELISA Kit	BSKM60804	96T
Mouse IL-1β ELISA Kit	bsk12015	96T
Mouse IL-1α ELISA Kit	bsk12025	96T
Mouse IL-2 ELISA Kit	bsk12016	96T
Mouse IL20 ELISA Kit	BSKM60806	96T
Mouse IL-21 ELISA Kit	bsk12026	96T
Mouse IL-22 ELISA Kit	bsk12034	96T
Mouse IL22Ra2 ELISA Kit	BSKM60807	96T
Mouse IL-23 ELISA Kit	bsk12029	96T
Mouse IL24 ELISA Kit	BSKM60808	96T
Mouse IL-25/IL-17E ELISA Kit	bsk12027	96T
Mouse IL-27 ELISA Kit	bsk12035	96T
Mouse IL28A ELISA Kit	BSKM60596	96T
Mouse IL2Ra ELISA Kit	BSKM60598	96T
Mouse IL3 ELISA Kit	BSKM60353	96T
Mouse IL31 ELISA Kit	BSKM60600	96T
Mouse IL-33 ELISA Kit	bsk12009	96T
Mouse IL34 ELISA Kit	BSKM60810	96T

产品名称	产品编号	产品规格
 Mouse		
Mouse IL-4 ELISA Kit	bsk12003	96T
Mouse IL-5 ELISA Kit	bsk12019	96T
Mouse IL-6 ELISA Kit	bsk12004	96T
Mouse IL6R ELISA Kit	BSKM60610	96T
Mouse IL-7 ELISA Kit	bsk12032	96T
Mouse IL-8/CXCL8 ELISA Kit	bsk12005	96T
Mouse IL8Ra ELISA Kit	BSKM60357	96T
Mouse IL8Rb ELISA Kit	BSKM60817	96T
Mouse IL-9 ELISA Kit	bsk12022	96T
Mouse Insulin (INS) ELISA Kit	BSKM61518	96T
Mouse Insulin Receptor Substrate 1 (IRS1) ELISA Kit	BSKM64861	96T
Mouse Interferon Gamma (IFNg) ELISA Kit	BSKM60211	96T
Mouse Interleukin 1 Beta (IL1b) ELISA Kit	BSKM60085	96T
Mouse Interleukin 10 (IL10) ELISA Kit	BSKM60084	96T
Mouse Interleukin 6 (IL6) ELISA Kit	BSKM60086	96T
Mouse IRF1 ELISA Kit	BSKM60617	96T
Mouse IRF3 ELISA Kit	BSKM60618	96T
Mouse IRF5 ELISA Kit	BSKM60619	96T
Mouse ITaC ELISA Kit	BSKM60212	96T
Mouse KLK1 ELISA Kit	BSKM60498	96T
Mouse Klotho Beta (KLb) ELISA Kit	BSKM64874	96T
Mouse Leptin ELISA Kit	bsk12012	96T
Mouse LIX ELISA Kit	bsk12078	96T
Mouse LOX1 ELISA Kit	BSKM60097	96T
Mouse LTb ELISA Kit	BSKM60620	96T
Mouse LTBP1 ELISA Kit	BSKM60820	96T
Mouse Luteinizing Hormone (LH) ELISA Kit	BSKM61699	96T
Mouse LY75 ELISA Kit	BSKM60621	96T
Mouse Macrophage Migration Inhibitory Factor (MIF) ELISA Kit	BSKM61547	96T
Mouse MCP-1 ELISA Kit	bsk12011	96T
Mouse MCP-2/CCL8 ELISA Kit	bsk12064	96T
Mouse MCP5 ELISA Kit	BSKM60158	96T
Mouse MDC ELISA Kit	BSKM60173	96T
Mouse MFG8 ELISA Kit	BSKM60041	96T
Mouse MIP-1 alpha/CCL3 ELISA Kit	bsk12042	96T
Mouse MIP-1 beta/CCL4 ELISA Kit	bsk12062	96T
Mouse MIP-2/CXCL2 ELISA Kit	bsk12024	96T
Mouse MIP3a ELISA Kit	BSKM60174	96T
Mouse MMP-3 ELISA Kit	bsk12068	96T
Mouse MMP-9 ELISA Kit	bsk12013	96T
Mouse MUC5B ELISA Kit	BSKM60499	96T
Mouse Mucin 2 (MUC2) ELISA Kit	BSKM61563	96T
Mouse Myeloid Differentiation Factor 88 (MyD88) ELISA Kit	BSKM62769	96T
Mouse Myeloperoxidase ELISA Kit	bsk12072	96T
Mouse Myogenin (MYOG) ELISA Kit	BSKM62771	96T

产品名称	产品编号	产品规格
 Mouse		
Mouse Myosin Heavy Chain 1, Skeletal Muscle, Adult (MYH1) ELISA Kit	BSKM61565	96T
Mouse Myosin Heavy Chain 2, Skeletal Muscle, Adult (MYH2) ELISA Kit	BSKM61566	96T
Mouse Myosin Heavy Chain 4, Skeletal Muscle (MYH4) ELISA Kit	BSKM64932	96T
Mouse NACHT, LRR and PYD Domains-Containing Protein 3(NLRP3) ELISA Kit	BSKM65170	96T
Mouse Neuropeptide Y (NPY) ELISA Kit	BSKM61578	96T
Mouse NFkB ELISA Kit	BSKM60501	96T
Mouse NGAL/Lipocalin-2 ELISA Kit	bsk12054	96T
Mouse NGF ELISA Kit	BSKM60213	96T
Mouse Nicotinamide-N-Methyltransferase (NNMT) ELISA Kit	BSKM64947	96T
Mouse Noggin (NOG) ELISA Kit	BSKM64949	96T
Mouse NT3 ELISA Kit	BSKM60502	96T
Mouse N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) ELISA Kit	BSKM61583	96T
Mouse OPN ELISA Kit	BSKM60093	96T
Mouse Osteocalcin (OC) ELISA Kit	BSKM61585	96T
Mouse Oxidized Low Density Lipoprotein (OxLDL) ELISA Kit	BSKM61591	96T
Mouse PD-1 ELISA Kit	bsk12070	96T
Mouse PDGFA ELISA Kit	BSKM60359	96T
Mouse PDGFB ELISA Kit	BSKM60822	96T
Mouse PDGFC ELISA Kit	BSKM60823	96T
Mouse PDGFRL ELISA Kit	BSKM60824	96T
Mouse PF4 ELISA Kit	BSKM60177	96T
Mouse Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide (PIK3C2a) ELISA Kit	BSKM61610	96T
Mouse PLGF ELISA Kit	BSKM60014	96T
Mouse PLGF-2 ELISA Kit	bsk12041	96T
Mouse Progesterone Receptor (PGR) ELISA Kit	BSKM62833	96T
Mouse Prolactin ELISA Kit	bsk12079	96T
Mouse Protein Kinase B Alpha (PKBa) ELISA Kit	BSKM64988	96T
Mouse Protein Kinase B Beta (PKBb) ELISA Kit	BSKM62841	96T
Mouse PTEN ELISA Kit	BSKM60503	96T
Mouse PTN ELISA Kit	BSKM60627	96T
Mouse RBP4 ELISA Kit	BSKM60179	96T
Mouse RETN ELISA Kit	BSKM60094	96T
Mouse S100 Calcium Binding Protein A8 (S100A8) ELISA Kit	BSKM62892	96T
Mouse S100A10 ELISA Kit	BSKM60504	96T
Mouse S100A6 ELISA Kit	BSKM60505	96T
Mouse SAP ELISA Kit	BSKM60154	96T
Mouse SCF ELISA Kit	BSKM60102	96T
Mouse Secretory Leukocyte Peptidase Inhibitor (SLPI) ELISA Kit	BSKM62911	96T
Mouse SELL ELISA Kit	BSKM60098	96T
Mouse SELP ELISA Kit	BSKM60099	96T
Mouse sICAM-1 ELISA Kit	bsk12055	96T
Mouse SIRT2 ELISA Kit	BSKM61275	96T
Mouse SLC ELISA Kit	BSKM60628	96T
Mouse Somatostatin (SST) ELISA Kit	BSKM61650	96T
Mouse Substance P (SP) ELISA Kit	BSKM61702	96T
Mouse Suppressors Of Cytokine Signaling 1 (SOCS1) ELISA Kit	BSKM65073	96T

产品名称	产品编号	产品规格
 Mouse		
Mouse Suppressors Of Cytokine Signaling 3 (SOCS3) ELISA Kit	BSKM62915	96T
Mouse Synapsin I (SYN1) ELISA Kit	BSKM65095	96T
Mouse Synaptophysin (SYP) ELISA Kit	BSKM61652	96T
Mouse TGFa ELISA Kit	BSKM60360	96T
Mouse TGFb1 ELISA Kit	BSKM60038	96T
Mouse TGFb2 ELISA Kit	BSKM60361	96T
Mouse TGFb3 ELISA Kit	BSKM60634	96T
Mouse TGFb1 ELISA Kit	BSKM60362	96T
Mouse TGFbR1 ELISA Kit	BSKM60363	96T
Mouse TGFbR2 ELISA Kit	BSKM60828	96T
Mouse Thrombomodulin/BDCA-3 ELISA Kit	bsk12033	96T
Mouse Thrombopoietin ELISA Kit	bsk12082	96T
Mouse Thrombospondin 2 (THBS2) ELISA Kit	BSKM65102	96T
Mouse Thymus Activation Regulated Chemokine (TARC) ELISA Kit	BSKM62930	96T
Mouse Thyroid Stimulating Hormone (TSH) ELISA Kit	BSKM61681	96T
Mouse Tie2 ELISA Kit	BSKM60113	96T
Mouse TIMP1 ELISA Kit	BSKM60208	96T
Mouse TIMP2 ELISA Kit	BSKM60235	96T
Mouse TL1A/TNFSF15 ELISA Kit	bsk12077	96T
Mouse TNFb ELISA Kit	BSKM60364	96T
Mouse TNFRSF17 ELISA Kit	BSKM60365	96T
Mouse TNFRSF4 ELISA Kit	BSKM60638	96T
Mouse TNFRSF5 ELISA Kit	BSKM60036	96T
Mouse TNFSF12 ELISA Kit	BSKM60830	96T
Mouse TNFSF14 ELISA Kit	BSKM60366	96T
Mouse TNF-α ELISA Kit	bsk12002	96T
Mouse TRAIL ELISA Kit	BSKM60367	96T
Mouse TRANCE/RANKL/TNFSF11 ELISA Kit	bsk12048	96T
Mouse TREM1 ELISA Kit	BSKM60124	96T
Mouse TRF ELISA Kit	BSKM60045	96T
Mouse Triggering Receptor Expressed On Myeloid Cells 2 (TREM2) ELISA Kit	BSKM65125	96T
Mouse Troponin I Type 3, Cardiac (TNNI3) ELISA Kit	BSKM61671	96T
Mouse Troponin T Type 2, Cardiac (TNNT2) ELISA Kit	BSKM65112	96T
Mouse TSLP ELISA Kit	bsk12040	96T
Mouse Tubulin Beta 3 (TUBb3) ELISA Kit	BSKM66476	96T
Mouse Vasoactive Intestinal Peptide (VIP) ELISA Kit	BSKM61691	96T
Mouse VCAM-1 ELISA Kit	bsk12039	96T
Mouse VEGF ELISA Kit	bsk12018	96T
Mouse VEGFB ELISA Kit	BSKM60368	96T
Mouse VEGFC ELISA Kit	BSKM60369	96T
Mouse VEGFD ELISA Kit	BSKM60370	96T
Mouse VEGFR1 ELISA Kit	BSKM60023	96T
Mouse VEGFR2 ELISA Kit	BSKM60029	96T
Mouse VEGFR3 ELISA Kit	BSKM66501	96T
Mouse VGF ELISA Kit	BSKM60646	96T
Mouse Wingless Type MMTV Integration Site Family, Member 3A (WNT3A) ELISA Kit	BSKM66330	96T
Mouse Wingless Type MMTV Integration Site Family, Member 5A (WNT5A) ELISA Kit	BSKM66332	96T

产品名称	产品编号	产品规格
 Rat		
Rat a2M ELISA Kit	BSKR60709	96T
Rat ACE2 ELISA Kit	BSKR60832	96T
Rat ADIPOR1 ELISA Kit	BSKR60506	96T
Rat ADIPOR2 ELISA Kit	BSKR60507	96T
Rat ADP ELISA Kit	BSKR60509	96T
Rat AGT ELISA Kit	BSKR60512	96T
Rat AGTR1 ELISA Kit	BSKR60833	96T
Rat AGTR2 ELISA Kit	BSKR60511	96T
Rat ALB ELISA Kit	BSKR60001	96T
Rat ALCAM ELISA Kit	BSKR60513	96T
Rat Amyloid Beta Peptide 1-42 (Ab1-42) ELISA Kit	BSKR62363	96T
Rat AngII ELISA Kit	BSKR60514	96T
Rat Angiotensin I Converting Enzyme (ACE) ELISA Kit	BSKR62367	96T
Rat Anti-Mullerian Hormone (AMH) ELISA Kit	BSKR62392	96T
Rat APOA1 ELISA Kit	BSKR60711	96T
Rat Atrial Natriuretic Peptide (ANP) ELISA Kit	BSKR62400	96T
Rat B-Cell Activating Factor (BAFF) ELISA Kit	BSKR65208	96T
Rat Bcl2 ELISA Kit	BSKR60517	96T
Rat BDNF ELISA Kit	BSKR60236	96T
Rat BLC1 ELISA Kit	BSKR60834	96T
Rat BMP2 ELISA Kit	BSKR60521	96T
Rat BMP4 ELISA Kit	BSKR60524	96T
Rat BMPR1A ELISA Kit	BSKR60526	96T
Rat bTG ELISA Kit	BSKR60230	96T
Rat Calcitonin Gene Related Peptide (CGRP) ELISA Kit	BSKR62499	96T
Rat CASP7 ELISA Kit	BSKR60713	96T
Rat CASP8 ELISA Kit	BSKR60714	96T
Rat Catalase (CAT) ELISA Kit	BSKR65649	96T
Rat Cholecystokinin (CCK) ELISA Kit	BSKR62475	96T
Rat Cluster Of Differentiation 40 Ligand (CD40L) ELISA Kit	BSKR62487	96T
Rat C-Peptide ELISA Kit	BSKR60715	96T
Rat CRP ELISA Kit	bsk13010	96T
Rat CTGF ELISA Kit	BSKR60527	96T
Rat Cyclin Dependent Kinase Inhibitor 1A (CDKN1A) ELISA Kit	BSKR65651	96T
Rat Diamine Oxidase (DAO) ELISA Kit	BSKR62544	96T
Rat EDN1 ELISA Kit	BSKR60529	96T
Rat EGF ELISA Kit	bsk13013	96T
Rat EGFR ELISA Kit	BSKR60532	96T
Rat EGFR2 ELISA Kit	BSKR60835	96T
Rat EG-VEGF ELISA Kit	BSKR60533	96T
Rat Elastase 2, Neutrophil (ELA2) ELISA Kit	BSKR62559	96T
Rat Enolase, Neuron Specific (NSE) ELISA Kit	BSKR62793	96T
Rat FABP2 ELISA Kit	BSKR60716	96T
Rat FAS ELISA Kit	BSKR60717	96T

产品名称	产品编号	产品规格
 Rat		
Rat FGF1 ELISA Kit	BSKR60536	96T
Rat FGF10 ELISA Kit	BSKR60836	96T
Rat FGF2 ELISA Kit	BSKR60237	96T
Rat FGF23 ELISA Kit	BSKR60537	96T
Rat FGF4 ELISA Kit	BSKR60541	96T
Rat FGF6 ELISA Kit	BSKR60542	96T
Rat FGF7 ELISA Kit	BSKR60543	96T
Rat FGF9 ELISA Kit	BSKR60545	96T
Rat FGFR1 ELISA Kit	BSKR60837	96T
Rat Fibrinogen Alpha (FGa) ELISA Kit	BSKR65255	96T
Rat Fibrinogen Beta (FGb) ELISA Kit	BSKR65256	96T
Rat Fibrinogen Gamma (FGg) ELISA Kit	BSKR65726	96T
Rat Follicle Stimulating Hormone (FSH) ELISA Kit	BSKR60303	96T
Rat GAL1 ELISA Kit	BSKR60131	96T
Rat Gastrin (GT) ELISA Kit	BSKR65273	96T
Rat GCSF ELISA Kit	BSKR60238	96T
Rat GDF2 ELISA Kit	BSKR60838	96T
Rat GDNF ELISA Kit	BSKR60718	96T
Rat Glucoprotein 130 (gp130) ELISA Kit	BSKR62619	96T
Rat Glucuronidase Beta (GUSb) ELISA Kit	BSKR65274	96T
Rat Glutamate Dehydrogenase 1 (GDH) ELISA Kit	BSKR62606	96T
Rat GM-CSF ELISA Kit	bsk13021	96T
Rat GROb ELISA Kit	BSKR60109	96T
Rat GROg ELISA Kit	BSKR60172	96T
Rat HBEGF ELISA Kit	BSKR60239	96T
Rat Hecpudin (Hepc) ELISA Kit	BSKR65277	96T
Rat HGF ELISA Kit	BSKR60548	96T
Rat HSPA8 ELISA Kit	BSKR60719	96T
Rat HSPb6 ELISA Kit	BSKR60720	96T
Rat Hypoxia Inducible Factor 1 Alpha (HIF1a) ELISA Kit	BSKR62642	96T
Rat ICAM-1/CD54 ELISA Kit	bsk13022	96T
Rat IFNa ELISA Kit	BSKR60088	96T
Rat IFNb ELISA Kit	BSKR60553	96T
Rat IFNg ELISA Kit	BSKR60058	96T
Rat IFN-γ ELISA Kit	bsk13002	96T
Rat IGF1R ELISA Kit	BSKR60839	96T
Rat IGF2 ELISA Kit	BSKR60558	96T
Rat IGFBP1 ELISA Kit	BSKR60560	96T
Rat IGFBP2 ELISA Kit	BSKR60562	96T
Rat IGFBP3 ELISA Kit	BSKR60563	96T
Rat IGFBP4 ELISA Kit	BSKR60564	96T
Rat IGFBP6 ELISA Kit	BSKR60565	96T
Rat IGF-I ELISA Kit	bsk13023	96T
Rat IgG total ELISA Kit	bsk13018	96T

产品名称	产品编号	产品规格
Rat		
Rat IgM ELISA Kit	bsk13019	96T
Rat IL-10 ELISA Kit	bsk13014	96T
Rat IL10Ra ELISA Kit	BSKR60840	96T
Rat IL10Rb ELISA Kit	BSKR60841	96T
Rat IL11 ELISA Kit	BSKR60571	96T
Rat IL12A ELISA Kit	BSKR60573	96T
Rat IL13 ELISA Kit	BSKR60576	96T
Rat IL15 ELISA Kit	BSKR60578	96T
Rat IL16 ELISA Kit	BSKR60580	96T
Rat IL17 ELISA Kit	BSKR60583	96T
Rat IL-17A ELISA Kit	bsk13009	96T
Rat IL17B ELISA Kit	BSKR60842	96T
Rat IL17F ELISA Kit	BSKR60843	96T
Rat IL18 ELISA Kit	BSKR60240	96T
Rat IL19 ELISA Kit	BSKR60587	96T
Rat IL1RA ELISA Kit	BSKR60590	96T
Rat IL-1β ELISA Kit	bsk13006	96T
Rat IL-1α ELISA Kit	bsk13017	96T
Rat IL-2 ELISA Kit	bsk13001	96T
Rat IL21 ELISA Kit	BSKR60844	96T
Rat IL23 ELISA Kit	BSKR60593	96T
Rat IL25 ELISA Kit	BSKR60845	96T
Rat IL2Ra ELISA Kit	BSKR60846	96T
Rat IL3 ELISA Kit	BSKR60602	96T
Rat IL33 ELISA Kit	BSKR60847	96T
Rat IL-4 ELISA Kit	bsk13004	96T
Rat IL5 ELISA Kit	BSKR60606	96T
Rat IL-6 ELISA Kit	bsk13005	96T
Rat IL6R ELISA Kit	BSKR60848	96T
Rat IL7 ELISA Kit	BSKR60611	96T
Rat IL7R ELISA Kit	BSKR60849	96T
Rat IL8Ra ELISA Kit	BSKR60612	96T
Rat IL9 ELISA Kit	BSKR60615	96T
Rat Insulin (INS) ELISA Kit	BSKR62670	96T
Rat Interleukin 1 Beta (IL1b) ELISA Kit	BSKR60104	96T
Rat IP10 ELISA Kit	BSKR60616	96T
Rat IRAK4 ELISA Kit	BSKR60850	96T
Rat IRF4 ELISA Kit	BSKR60851	96T
Rat ITaC ELISA Kit	BSKR60242	96T
Rat KLK1 ELISA Kit	BSKR60721	96T
Rat Leptin (LEP) ELISA Kit	BSKR62700	96T
Rat LFA1a ELISA Kit	BSKR60852	96T
Rat LTb ELISA Kit	BSKR60853	96T
Rat Luteinizing Hormone (LH) ELISA Kit	BSKR63005	96T

产品名称	产品编号	产品规格
Rat		
Rat MCP-1 ELISA Kit	bsk13024	96T
Rat Mlg ELISA Kit	BSKR60854	96T
Rat MIP1a ELISA Kit	BSKR60304	96T
Rat MIP-2 ELISA Kit	bsk13015	96T
Rat Monocyte Chemotactic Protein 1 (MCP1) ELISA Kit	BSKR62723	96T
Rat Motilin (MTL) ELISA Kit	BSKR62763	96T
Rat Myeloperoxidase (MPO) ELISA Kit	BSKR62755	96T
Rat N-Acetyl Beta-D-Glucosaminidase (NAGase) ELISA Kit	BSKR62774	96T
Rat NAP3 ELISA Kit	BSKR60118	96T
Rat Nerve Growth Factor (NGF) ELISA Kit	BSKR60243	96T
Rat NGAL/Lipocalin-2 ELISA Kit	bsk13020	96T
Rat N-Terminal Pro Brain Natriuretic Peptide (NT-ProBNP) ELISA Kit	BSKR62800	96T
Rat Osteocalcin (OC) ELISA Kit	BSKR62804	96T
Rat Osteopontin (OPN) ELISA Kit	BSKR62810	96T
Rat Oxidized Low Density Lipoprotein (OxLDL) ELISA Kit	BSKR62813	96T
Rat PGA ELISA Kit	BSKR60722	96T
Rat Phospholipase A2, Lipoprotein Associated (LpPLA2) ELISA Kit	BSKR62707	96T
Rat PLGF ELISA Kit	BSKR60626	96T
Rat PTEN ELISA Kit	BSKR60723	96T
Rat PTN ELISA Kit	BSKR60855	96T
Rat Retinol Binding Protein 4, Plasma (RBP4) ELISA Kit	BSKR62881	96T
Rat RETN ELISA Kit	BSKR60724	96T
Rat S100A10 ELISA Kit	BSKR60725	96T
Rat S100A6 ELISA Kit	BSKR60726	96T
Rat S100A8 ELISA Kit	BSKR60727	96T
Rat SELP ELISA Kit	BSKR60728	96T
Rat Stem Cell Factor (SCF) ELISA Kit	BSKR62899	96T
Rat Stromal Cell Derived Factor 1 (SDF1) ELISA Kit	BSKR62902	96T
Rat Substance P (SP) ELISA Kit	BSKR63007	96T
Rat Suppressors Of Cytokine Signaling 1 (SOCS1) ELISA Kit	BSKR65923	96T
Rat Suppressors Of Cytokine Signaling 3 (SOCS3) ELISA Kit	BSKR65358	96T
Rat TGFa ELISA Kit	BSKR60629	96T
Rat TGFb1 ELISA Kit	BSKR60631	96T
Rat TGFb2 ELISA Kit	BSKR60633	96T
Rat TGFb3 ELISA Kit	BSKR60856	96T
Rat Thymus Activation Regulated Chemokine (TARC) ELISA Kit	BSKR65367	96T
Rat Thyroid Stimulating Hormone (TSH) ELISA Kit	BSKR62969	96T
Rat TIMP-1 ELISA Kit	bsk13012	96T
Rat TIMP2 ELISA Kit	BSKR60305	96T
Rat TNFb ELISA Kit	BSKR60636	96T
Rat TNFRSF1A ELISA Kit	BSKR60857	96T
Rat TNFRSF1B ELISA Kit	BSKR60858	96T
Rat TNFRSF5 ELISA Kit	BSKR60859	96T
Rat TNF-α ELISA Kit	bsk13003	96T

产品名称	产品编号	产品规格
------	------	------

Rat

Rat Total MMP-9 ELISA Kit	bsk13011	96T
Rat TRAIL ELISA Kit	BSKR60640	96T
Rat TRF ELISA Kit	BSKR60306	96T
Rat Troponin T Type 2, Cardiac (TNNT2) ELISA Kit	BSKR65945	96T
Rat TSLP ELISA Kit	BSKR60860	96T
Rat Vasoactive Intestinal Peptide (VIP) ELISA Kit	BSKR62990	96T
Rat VEGF ELISA Kit	bsk13007	96T
Rat VEGFA ELISA Kit	BSKR60244	96T
Rat VEGFB ELISA Kit	BSKR60644	96T
Rat VEGFC ELISA Kit	BSKR60729	96T
Rat VEGFD ELISA Kit	BSKR60645	96T
Rat VEGFR1 ELISA Kit	BSKR60861	96T
Rat VEGFR2 ELISA Kit	BSKR60862	96T
Rat VEGFR3 ELISA Kit	BSKR60863	96T
Rat VGF ELISA Kit	BSKR60864	96T
Rat β -NGF ELISA Kit	bsk13008	96T

Rabbit

Rabbit Interleukin 4 (IL4) ELISA Kit	BSKO60604	96T
Rabbit Interleukin 6 (IL6) ELISA Kit	BSKO60609	96T
Rabbit Matrix Metalloproteinase 13 (MMP13) ELISA Kit	BSKO62742	96T
Rabbit Tumor Necrosis Factor Alpha (TNFa) ELISA Kit	BSKO60108	96T

Chicken

Chicken Interferon Gamma (IFNg) ELISA Kit	BSKC60554	96T
Chicken Interleukin 1 Beta (IL1b) ELISA Kit	BSKC60589	96T
Chicken Interleukin 10 (IL10) ELISA Kit	BSKC60567	96T
Chicken Interleukin 12A (IL12A) ELISA Kit	BSKC60572	96T
Chicken Interleukin 12B (IL12B) ELISA Kit	BSKC60574	96T
Chicken Interleukin 13 (IL13) ELISA Kit	BSKC60575	96T
Chicken Interleukin 16 (IL16) ELISA Kit	BSKC60579	96T
Chicken Interleukin 18 (IL18) ELISA Kit	BSKC60584	96T
Chicken Interleukin 2 (IL-2) ELISA Kit	BSKC60597	96T
Chicken Interleukin 4 (IL4) ELISA Kit	BSKC60603	96T
Chicken Interleukin 6 (IL6) ELISA Kit	BSKC60608	96T
Chicken Interleukin 8 (IL8) ELISA Kit	BSKC60241	96T
Chicken Nuclear Factor Kappa B (NFkB) ELISA Kit	BSKC65317	96T
Chicken Transforming Growth Factor Beta 1 (TGFb1) ELISA Kit	BSKC60630	96T
Chicken Vascular Endothelial Growth Factor A (VEGFA) ELISA Kit	BSKC60642	96T

产品名称	产品编号	产品规格
------	------	------

Swine

Swine IFN- γ ELISA Kit	bsk14006	96T
Swine IL-2 ELISA Kit	bsk14003	96T
Swine IL-4 ELISA Kit	bsk14004	96T
Swine IL-6 ELISA Kit	bsk14001	96T
Swine IL-8 ELISA Kit	bsk14008	96T
Swine IL-10 ELISA Kit	bsk14005	96T
Swine IL-12/IL-23p40 ELISA Kit	bsk14009	96T
Swine IL-1 β ELISA Kit	bsk14002	96T
Swine TNF- α ELISA Kit	bsk14007	96T

Porcine

Pig IL8 ELISA Kit	BSKP60481	96T
Pig TNFa ELISA Kit	BSKP60483	96T
Porcine Epidermal Growth Factor (EGF) ELISA Kit	BSKP61716	96T
Porcine Insulin Like Growth Factor 1 (IGF1) ELISA Kit	BSKP61733	96T
Porcine Interleukin 17 (IL17) ELISA Kit	BSKP61756	96T
Porcine Interleukin 2 (IL2) ELISA Kit	BSKP61766	96T
Porcine Interleukin 4 (IL4) ELISA Kit	BSKP61768	96T
Porcine Interleukin 6 (IL6) ELISA Kit	BSKP61772	96T
Porcine Troponin I Type 3, Cardiac (TNNI3) ELISA Kit	BSKP65578	96T
Porcine Troponin T Type 2, Cardiac (TNNT2) ELISA Kit	BSKP66463	96T

Bovine

Bovine Beta-Lactoglobulin (bLg) ELISA Kit	BSKB62327	96T
Bovine Fibronectin (FN) ELISA Kit	BSKB60162	96T
Bovine Interferon Alpha (IFNa) ELISA Kit	BSKB61727	96T
Bovine Interferon Beta (IFNb) ELISA Kit	BSKB61728	96T
Bovine Interferon Gamma Induced Protein 10kDa (IP10) ELISA Kit	BSKB61773	96T
Bovine Interleukin 1 Beta (IL1b) ELISA Kit	BSKB61761	96T
Bovine Interleukin 6 (IL6) ELISA Kit	BSKB61771	96T
Bovine Interleukin 8 (IL8) ELISA Kit	BSKB60223	96T
Bovine Interleukin 18 (IL18) ELISA Kit	BSKB61757	96T
Bovine Pepsinogen A (PGA) ELISA Kit	BSKB65531	96T
Bovine Regulated On Activation In Normal T-Cell Expressed And Secreted (RANTES) ELISA Kit	BSKB65543	96T
Bovine Tumor Necrosis Factor Alpha (TNFa) ELISA Kit	BSKB60482	96T

产品名称	产品编号	产品规格
------	------	------

Guinea pig

Guinea pig C Reactive Protein (CRP) ELISA Kit	BSKI64703	96T
Guinea pig Interferon Gamma (IFN γ) ELISA Kit	BSKI60790	96T
Guinea pig Interleukin 1 Beta (IL1 β) ELISA Kit	BSKI60187	96T
Guinea pig Interleukin 2 (IL2) ELISA Kit	BSKI64852	96T
Guinea pig Interleukin 4 (IL4) ELISA Kit	BSKI60813	96T

Monkey

Monkey Insulin (INS) ELISA Kit	BSKA66401	96T
Monkey Interferon Gamma (IFN γ) ELISA Kit	BSKA64557	96T
Monkey Interleukin 1 Beta (IL1 β) ELISA Kit	BSKA66400	96T
Monkey Interleukin 6 (IL6) ELISA Kit	BSKA64563	96T

Canine

Canine Interferon Gamma (IFN γ) ELISA Kit	BSKD60789	96T
Canine Interleukin 1 Beta (IL1 β) ELISA Kit	BSKD60186	96T
Canine Interleukin 4 (IL-4) ELISA Kit	BSKD60812	96T
Canine Osteocalcin (OC) ELISA Kit	BSKD64963	96T
Canine Tumor Necrosis Factor Alpha (TNF α) ELISA Kit	BSKD60193	96T

SARS-CoV-2

SARS-CoV-2 (2019-nCoV) IgM Antibody ELISA Kit	BSKV0002	96T
SARS-CoV-2 (2019-nCoV) Neutralization Antibody Screening ELISA Kit	BSKV0004	96T
SARS-CoV-2 (2019-nCoV) Nucleoprotein ELISA Kit	BSKV0001	96T
SARS-CoV-2 (2019-nCoV) Total Antibody Detection ELISA Kit	BSKV0005	96T
SARS-CoV-2 (2019-nCoV) Total Antibody ELISA Kit	BSKV0003	96T

产品名称	产品编号	产品规格
------	------	------

Other

Aflatoxin B1 ELISA Kit	bsk60001	96T
Aflatoxin B1 ELISA Kit	bsk60002	96T
Aflatoxin B1 ELISA Kit	bsk60003	96T
Aflatoxin B1 ELISA Kit	bsk60004	96T
Brucellosis IgM ELISA Kit	bsk31013	96T
Footand Mouth Disease Virus ELISA Kit	bsk31014	96T
General 17-OHP ELISA Kit	BSKL60730	96T
General 5-HT ELISA Kit	BSKL60731	96T
General 8-OHdG ELISA Kit	BSKL60732	96T
General Acetylcholine (ACH) ELISA Kit	BSKL66334	96T
General Adenosine Diphosphate (ADP) ELISA Kit	BSKL61342	96T
General Adenosine Triphosphate (ATP) ELISA Kit	BSKL66336	96T
General Advanced Glycation End Product (AGE) ELISA Kit	BSKL64554	96T
General Aldosterone (ALD) ELISA Kit	BSKL66335	96T
General Arachidonic Acid (AA) ELISA Kit	BSKL66385	96T
General BK ELISA Kit	BSKL66470	96T
General cAMP ELISA Kit	BSKL60246	96T
General CAR ELISA Kit	BSKL61344	96T
General Cor ELISA Kit	BSKL60737	96T
General Cort ELISA Kit	BSKL60736	96T
General Creatinine ELISA Kit	BSKL60738	96T
General Cyclic Guanosine Monophosphate (cGMP) ELISA Kit	BSKL66337	96T
General DDAVP ELISA Kit	BSKL61345	96T
General DHEA ELISA Kit	BSKL60739	96T
General DHVD3 ELISA Kit	BSKL60741	96T
General Dopamine (DA) ELISA Kit	BSKL66339	96T
General E1 ELISA Kit	BSKL61346	96T
General E1S ELISA Kit	BSKL60742	96T
General E2 ELISA Kit	BSKL60308	96T
General E3 ELISA Kit	BSKL60743	96T
General EPI ELISA Kit	BSKL60744	96T
General ET ELISA Kit	BSKL60745	96T
General Gamma-Aminobutyric Acid (gABA) ELISA Kit	BSKL62329	96T
General GC ELISA Kit	BSKL60748	96T
General HA ELISA Kit	BSKL60749	96T
General Histamine (HIS) ELISA Kit	BSKL66342	96T
General HVD3 ELISA Kit	BSKL60750	96T
General Hydroxyproline (Hyp) ELISA Kit	BSKL66343	96T
General Inositol Triphosphate (IP3) ELISA Kit	BSKL62355	96T
General Leukotriene B4 (LTB4) ELISA Kit	BSKL66344	96T
General LPS ELISA Kit	BSKL61347	96T
General LXA4 ELISA Kit	BSKL60766	96T
General Melatonin (MT) ELISA Kit	BSKL66345	96T
General NE ELISA Kit	BSKL60751	96T

产品名称	产品编号	产品规格
Other		
General OT ELISA Kit	BSKL61348	96T
General Pg ELISA Kit	BSKL60755	96T
General PGH2 ELISA Kit	BSKL60756	96T
General PG12 ELISA Kit	BSKL60757	96T
General Platelet Activating Factor (PAF) ELISA Kit	BSKL66346	96T
General Prostaglandin E2 (PGE2) ELISA Kit	BSKL60753	96T
General Prostaglandin F2 Alpha (PGF2a) ELISA Kit	BSKL60754	96T
General PTD ELISA Kit	BSKL60758	96T
General SCT ELISA Kit	BSKL60759	96T
General Testosterone (T) ELISA Kit	BSKL60760	96T
General THB ELISA Kit	BSKL60761	96T
General Thyroxine (T4) ELISA Kit	BSKL60310	96T
General Triiodothyronine (T3) ELISA Kit	BSKL60309	96T
General TXA2 ELISA Kit	BSKL61349	96T
General TXB2 ELISA Kit	BSKL60762	96T
General VB12 ELISA Kit	BSKL60763	96T
General VB9 ELISA Kit	BSKL60764	96T
General VD2 ELISA Kit	BSKL66348	96T
General VK1 ELISA Kit	BSKL60765	96T
Ochratoxins A ELISA Kit	bsk60005	96T

* 产品信息均以官网为准: www.biosschina.com

试剂盒组成:

- 抗体预包被可拆 96 孔板
- 显色底物 (TMB), 反应终止液
- 酶标试剂
- 浓缩洗涤液, 样品稀释液
- 冻干标准品
- 封板胶纸, 产品说明书

* 不同检测指标试剂盒其组分可能有差距, 请按实际说明书中的组分列表为准。



07 / ELISA 配套试剂

ELISA 配套试剂列表

编号	名称	规格
C04-01001	ELISA 包被液 (1×)	100mL/500mL
C04-01002	ELISA 封闭液	100mL/500mL/1L
C2075	ELISA 封闭液 (1% 明胶)	100mL
C04-01003	ELISA 终止液	30mL/100mL
C04-01004	ELISA 洗涤液 (10X)	50mL/100mL
C04-01005	ELISA 酶标抗体稀释液	30mL/50mL/100mL
C04-02001	ELISA 酶标板	1 块 /96T
C04-02002	封板膜	1 张
C04-03001	单组分超灵敏 TMB 显色液	20mL/50mL/100mL
C04-03002	双组分 TMB 显色液	10mL/30mL/100mL

ELISA 配套试剂优势



品种完善



产品使用, 方便快捷



价格便宜质量好,
性价比高



售后完善,
轻松换货无烦恼

08 / ELISA 试剂盒定制

您可以提供特异性抗体，高纯度抗原，北京博奥森可根据您的个性化需要开发特定 ELISA 试剂盒。可采用 ELISA 夹心法、竞争法（测定抗原）、间接法（测定抗体）或竞争法（测定小分子化合物）定量检测蛋白的含量。

技术优势



服务内容

服务项目（编号：SV8000）	内容	价格	周期	交付结果
客户提供抗原、抗体原料，博奥森负责 ELISA KIT 组装 (5kits)	优化抗原、抗体配对条件、优化试剂盒、提供 5 个完整的试剂盒、标准曲线报告	询价	2-3 月	完整的、真实的、全套的检测数据报告及成型试剂盒。
客户提供抗原、抗体信息博奥森负责 ELISA KIT 开发	组合蛋白表达、抗体制备、蛋白标记，抗体对筛选，提供批内、批间差报告	询价	协商	

09 / ELISA 常见问题解答指南

可能原因	提供建议
------	------

问题一：高背景或阴性对照值偏高

阴性对照被阳性对照或样品污染	洗涤时，勿将洗液溢出孔外；若操作失误，请重新洗涤
洗涤不充分	确保在洗涤时微孔均充满液体，且将残余液体从孔中完全移除
酶结合物使用浓度过高	请检查酶结合物是否按说明书规定稀释
孵育温度过高	检查孵育箱的温度是否正确、稳定
孵育时间过长	请按照说明书规定时间孵育，误差要控制在 $\pm 1.5\text{min}$
底物在使用前曝光	使用前应避光保存
读板前停留时间过长	加终止液后 3 分钟内读数

问题二：吸光度（OD 值）偏低

试剂未平衡至室温	实验前将试剂盒从冰箱中取出并平衡至室温（提前 30min 即可）
孵育时间过短	请使用计时器以确保准确孵育时间
试剂在使用前未混匀	使用前将试剂混匀（3-5min）
酶标试剂使用浓度偏低	请按照说明书稀释至建议使用浓度
显色时间不足	适当延长显色时间

问题三：边缘效应

蒸发	整个实验过程必须使用封板胶纸密封微孔板
温度不均匀	校准孵育箱，勿叠放反应板

问题四：标准曲线不佳

倍比稀释标准品时未混匀	稀释标准品的每一步均需混匀（3-5min）
过早稀释	标准品在实验使用事前 10min 稀释即可
加入的体积不正确	使用校准的移液器，将标准品分配到各个微孔中

问题五：标准曲线良好，但样本无检测信号

样本中无检测物或检测物含量极低	设置内参，重新实验
样本中其他物质影响 / 掩盖检测	适当稀释以降低其他物质干扰，或作系列稀释，检测回收率
样本待检物浓度过高，Hook 效应致假阴性	适当倍数稀释样本，检测物浓度尽量在试剂盒检测范围内

问题六：重复性较差

微孔中有微小气泡	用注射器针头挑破气泡
试剂未充分混匀	确保充分混匀试剂（3-5min）
样本中有杂质或沉淀物	使用前离心
使用已用过的封板胶纸	每次洗板之后须使用新的封板胶纸封住反应孔



博奥森生物公众号



博奥森小红书账号



博奥森B站账号

北京博奥森生物技术有限公司

BEIJING BIOSYNTHESIS BIOTECHNOLOGY CO.,LTD.

🌐 www.biosschina.com ☎ 400-901-9800

🏠 北京市经济技术开发区景盛南四街联东U谷西区2号院67号楼